

LogAndStream for Shimmer3
Firmware User Manual
Rev 0.1a



Legal Notices and Disclaimer

Redistribution IS permitted provided that the following conditions are met:

Redistributions must retain the copyright notice, and the following disclaimer. Redistributions in electronic form must reproduce the above copyright notice, this list of conditions and the following disclaimer in the documentation and/or other materials provided with the document.

Neither the name of Shimmer Research, or Realtime Technologies Ltd. nor the names of its contributors may be used to endorse or promote products derived from this document without specific prior written permission.

THIS DOCUMENT IS PROVIDED BY THE COPYRIGHT HOLDERS AND CONTRIBUTORS "AS IS" AND ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE ARE DISCLAIMED. IN NO EVENT SHALL THE COPYRIGHT OWNER OR CONTRIBUTORS BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR CONSEQUENTIAL DAMAGES (INCLUDING, BUT NOT LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS OR SERVICES; LOSS OF USE, DATA, OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER CAUSED AND ON ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, STRICT LIABILITY, OR TORT (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF THE USE OF THIS DOCUMENT, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGE.



Table of Contents

| 1. | Intr | oduction | 3 |
|----|-----------------|---|------|
| | 1.1. | Scope of this User Manual | 3 |
| | 1.2. | Pre-Requisites | 3 |
| | 1.3. | Installation | 3 |
| 2. | Firn | nware Features | 4 |
| | 2.1. | General operation | 5 |
| | 2.2. | LED indicators | 6 |
| | 2.3. | Start/Stop Logging | 8 |
| 3. | | figuring the Shimmer3 | |
| | 3.1. | Configuration and calibration loading priority | |
| | 3.2. | Configuration file | |
| | 3.3. | Enabling sensors | |
| | 3.4. | Configuring sensors | |
| | 3.5. | Configuring an experiment | 14 |
| | 3.6. | Infomem contents | 15 |
| 4. | Blu | etooth Streaming | 18 |
| | 4.1. | Set Commands | 18 |
| | 4.2. | Get Commands | 18 |
| | 4.3. | Action Commands | 20 |
| | 4.4. | Streaming | 20 |
| | 4.5. | Bluetooth Latency | 21 |
| 5. | | ding the SD card Data | |
| | 5.1. | SD card directory structure | |
| | 5.2. | Parsing the raw data file | 22 |
| 6. | App 6.1. | Dendices Example of the <i>Shimmer3</i> configuration file (sdlog.cfg) | |
| | 6.2. | SD data file configuration header - sensor calibration bytes | 28 |
| | 6.3. | LogAndStream source code header file (Shimmer.h) | 31 |
| | 6.4 | Troubleshoot | . 41 |



1. Introduction

This document is an accompaniment to the *LogAndStream Firmware v0.1.0* image and later for *Shimmer3*. No previous development experience is required. **Note:** *LogAndStream Firmware v0.1.0* is a Beta release.

LogAndStream firmware is a complete data recording solution which merges features from our previous Shimmer3 firmware releases. The LogAndStream firmware facilitates logging of data from a Shimmer3 to the on-board SD card while also providing the ability to simultaneously stream data wireless connection to a Bluetooth-enabled PC. The firmware allows for full user configuration of the Shimmer3 using a configuration file stored on the SD card, or over Bluetooth using the provided ShimmerCapture PC application.

It is recommended that the *LogAndStream* firmware be used in conjunction with the *ShimmerCapture* and/or *ShimmerLog* PC applications, which are available for download from the Shimmer website (www.shimmersensing.com).

1.1. Scope of this User Manual

The purpose of this *User Manual* is to guide the user through the features of the *LogAndStream* firmware and to provide the required instructions to configure the data logging options, configure the data streaming options and to parse the recorded and received data. The *User Manual* does not provide an extensive explanation of the source code for the firmware.

1.2. Pre-Requisites

The *LogAndStream* firmware can be used with a *Shimmer3* device that has a connected microSD card. For *Shimmer3* compatibility, the chosen microSD card must be less than 32GB in capacity, not be an XC (eXtended Capacity) card and it must support 1-bit SPI mode.

For a wireless connection to the *Shimmer3*, a Bluetooth enabled PC and the *ShimmerCapture* is required to interface with *Shimmer3* and receive the streamed data. The *ShimmerCapture* application allows a user to fully configure the *Shimmer3* whilst providing the ability to save streaming data locally to the PC in real-time. Concurrently, a *Shimmer Dock* is required to allow access the microSD card on the Shimmer from a PC. The *Shimmer Dock* provides an alternative method for configuring the *Shimmer3* logging preferences (using the *ShimmerLog* application) whilst also facilitating data transfer (either RAW or calibrated) from the microSD card. Please note that legacy (black) Shimmer docks are not suitable for this purpose and a newer (white) dock is required.

For more information on the *Shimmer3* platform, please refer to the *Shimmer3 User Manual* which is available for download at www.shimmersensing.com.

1.3. Installation

Install the *LogAndStream Firmware* v0.1.0 image (LogAndStream_v0.1.0_shimmer3.txt) onto a *Shimmer3* device, using the *Shimmer3 Bootstrap Loader* (*Shimmer3 BSL*) application. The Shimmer3BSL.exe application is available via download from the Members area at www.shimmersensing.com.



2. Firmware Features

As stated previously, the *LogAndStream* firmware provides a seamless integration of the existing *Shimmer3* firmware releases of *SDLog* (*i.e.*, logging *Shimmer3* data to an on-board microSD card) and *BtStream* (*i.e.*, streaming *Shimmer3* data back to "host" device). As such, many users will already be familiar with the operation of the firmware and many its features. The *LogAndStream* firmware has been designed such that a user can still utilise the *Shimmer3* in a purely *BtStream* or *SDLog* mode if they desire. What sets the *LogAndStream* firmware apart from previous firmware releases is its ability to perform both *SDLog* and *BtStream* operations simultaneously.

A *Shimmer3* programmed with *LogAndStream* firmware can be in one of five states: *Idle, BT Connected, BT Streaming, BT Streaming + SD Logging* or *SD Logging* - as shown in Figure 4-1. The active states form two operational branches from which the user can choose to operate the *Shimmer3* device - one initiated by a Bluetooth connection (blue shaded area in Figure 4-1) and the other based on SD Logging operation (orange shaded area in Figure 4-1).




Figure 2-1 LogAndStream firmware operational hierarchy.

The key feature of the LogAndStream firmware can be viewed at the intersection of the two active branches (green shaded area in Figure 4-1). From a "host" application, the user can select to stream Shimmer3 sensor data over a Bluetooth connection at the same time as saving it locally to the on-board microSD card. An instance where this is of key benefit is the case where, in the process of data collection, a Shimmer3 might travel out of range from the "host" device. In this circumstance, the Shimmer3 will continue to save data to the microSD card such that there will be no break in the



data recording. Similarly, during normal SD logging operation, continuous Bluetooth streaming is not always desirable or required (e.g., during long data recordings where a lower battery consumption is key). By combining Bluetooth support with SD logging capability, the new LogAndStream firmware now facilitates the ability for the user to check the recording activity of a Shimmer3 at any stage during SD logging operation - without interrupting the data recording.

2.1. General operation

State 1 - Idle

When a *Shimmer3*, programmed with the *LogAndStream* firmware, is first powered on or reset, it enters the *Idle* state. It will remain in the *Idle* state until either a connection is made over the Bluetooth link (*i.e.*, by opening a serial connection), or the user instigates a single button press which will initiate *SD logging* mode. It should be noted that this operational case is valid when "Undock Start" is disabled in the *Shimmer3* configuration - see Section 3.5. When "Undock Start" is enabled, the *Shimmer3* will progress directly from *Idle* mode into *SD Logging* mode after a power-on or reset condition.

Note: To use the Bluetooth features of the *LogAndStream* firmware, the *Shimmer3* device must first be paired with a "host" device (*e.g.*, a PC), as outlined in the *Shimmer3 User Manual*.

- Path (a): To enter *BT Connected* state, a Bluetooth connection must be initiated by a "host" side application (*e.g.*, *ShimmerCapture*). A subsequent disconnection of the Bluetooth connection will send the *Shimmer3* back to the *Idle* state.
- Path (c): If "Undock Start" is disabled, a short button press on the Shimmer3 will initiate SD Logging state. Similarly, another short button press will terminate SD Logging mode and return the Shimmer3 to Idle state.

State 2 - BT Connected

In the *BT Connected* state, the *Shimmer3* can process various commands to configure its sensors and sampling parameters, set calibration parameters, send configuration settings back to the "host" (PC, mobile or other) and start sampling. When the *Shimmer3* is in the *BT Connected*, *BT Streaming* or *BT Streaming* + SD Logging states, there can be active communication between the *Shimmer3* and the host over the Bluetooth serial connection. Packets of bytes are sent in both directions and these can consist of commands, responses or data.

Path (b): The host application (e.g., ShimmerCapture) allows the user to fully configure the Shimmer3 and choose between either BT Streaming or BT Streaming + SD Logging recording modes. Alternatively, the short button press action can be used to start or stop BT Streaming + SD Logging mode from the initial BT Connected branch shown in Figure 2-1.



State 3 - BT Streaming

When a command to start *BT Streaming* is received, the *Shimmer3* enters the *BT Streaming* state and starts sampling data from its sensors and sending that data over the Bluetooth link. This state will continue until a command to stop logging is received, whereupon the *Shimmer3* returns to the *BT Connected* state. When in the *BT Streaming* state, a short button press can also be used to halt streaming and return the *Shimmer3* to *BT Connected* state. Closing the serial connection will put the *Shimmer3* back into the *Idle* state. Similarly, if the *Shimmer3* goes out-of-range from the "host" device for a certain period ($\approx 150 \text{ s}$), the *Shimmer3* will also return to *Idle* state.

State 4 - BT Streaming + SD Logging

As described previously, the *BT Streaming + SD Logging* state can be initiated either by using the "host" application (*e.g.*, *ShimmerCapture*) or using a short button press from the initial *BT Connected* state. When the *Shimmer3* enters the *BT Streaming+ SD Logging* state, the *Shimmer3* will start sampling data from its sensors and send this data over the Bluetooth link as well as store it locally in raw file format to the microSD card (see Section 5).

The *Shimmer3* will continue in this mode until either one of two conditions. Firstly, if a stop logging command is received from the "host" or, as mentioned previously, the user instigates a short button press, the *Shimmer3* will return to the *BT Connected* state. Secondly, in the event of either a Bluetooth disconnection, the serial port being closed or the *Shimmer3* travels out of range for a certain period ($\approx 150 \text{ s}$), the *Shimmer3* will automatically enter the *SD Logging* state (as illustrated by "Path (d)" in Figure 4-1).

State 5 - SD Logging

In the *SD Logging* state, the *Shimmer3* will continuously read data from its sensors and store this data locally to a raw data file in its on-board microSD card (for more information see Section 5). The *Shimmer3* will remain open to a Bluetooth connection by a host device whereupon, the Shimmer3 which automatically switch to the *BT Streaming + SD Logging* state without interruption to the SD logging data stream.

2.2. LED indicators

The *Shimmer3* has five LEDs in two locations: the lower LED location (location A) which contains the green (b), yellow¹ and red LEDs; and the upper LED location (location B) which contains the green (a) and blue LEDs.

¹ Note that what is referred to as the yellow LED may appear orange to some users.




Figure 2-2 Shimmer3 in enclosure.

The two upper LEDs are used to indicate state of operation for the *LogAndStream* firmware. The status of each depends on whether the *Shimmer3* is in a docked or undocked state, as shown in Table 2-1 below.

| | UNDOCKED | | DOCKED | |
|-------------------|--------------------|--------------------|--------------------|--------------------|
| LED | Green (b) | Blue | Green (b) | Blue |
| Power OFF | OFF | OFF | OFF | OFF |
| Idle mode | OFF | 0.1 s ON/2 s OFF | N/A | 0.1 s ON/2 s OFF |
| Connected | OFF | ON | N/A | ON |
| Streaming only | OFF | 1 s ON/1 s OFF | N/A | 1 s ON/1 s OFF |
| Logging only | 1 s ON/1 s OFF | OFF | N/A | N/A |
| Streaming+Logging | 1 s each | | N/A | N/A |
| | **** | ***** | | |
| Configuring OFF | | 0.1 s ON/0.1 s OFF | N/A | 0.1 s ON/0.1 s OFF |
| Error* | 0.1 s ON/0.1 s OFF |

Table 2-1 - LED indicators specifc to LogAndStream firmware (undocked and docked). * In Error mode, the Blue and Green (b) LEDs alternate.

Note: The Shimmer unit should never be placed in the dock while the operation LEDs indicate that it is configuring as this may cause a file-system error. Once configuration has begun, you must power off or reset the Shimmer unit before docking.

Note: It is not recommend to place the *Shimmer3* unit in the dock while it is logging data as this can cause SD card access problems. Once logging has begun, you must power off or reset the *Shimmer3* unit before docking.

The lower three LEDs are exclusively used across all *Shimmer3* firmware as a battery charge indicator and will show the colour corresponding to the charge status. As above, the status of each LED depends on whether the *Shimmer3* is in a docked or undocked state. If the *Shimmer3* is on a dock or



multi-charger and is powered on, the battery charge status indicator will show a solid LED. If the *Shimmer3* is undocked, the corresponding LED will flash, as described below in Table 2-2.

| | | DOCKED | | | | |
|---------------|------------------|------------------|------------------|-----------|--------|-----|
| LED | Green (a) | Yellow | Red | Green (a) | Yellow | Red |
| Power OFF | OFF | OFF | OFF | OFF | OFF | OFF |
| Low Charge | OFF | OFF | 0.1 s ON/5 s OFF | OFF | OFF | ON |
| Medium Charge | OFF | 0.1 s ON/5 s OFF | OFF | OFF | ON | OFF |
| Full Charge | 0.1 s ON/5 s OFF | OFF | OFF | ON | OFF | OFF |

Table 2-2 - Power status LED indicators for Shimmer3 firmware (undocked and docked).

2.3. Start/Stop Logging

Dock (automatic start)

This is the default way for the *Shimmer3* to start and stop logging data to the SD card and is enabled if the "Undock Start" option is enabled from ShimmerLog or ShimmerCapture (*i.e.*, configuration file line entry: "userbutton=0"). With this configuration, if the Shimmer is powered on or reset while it is **not** on the dock, the following process will begin immediately:

- 1. The Shimmer unit will go into standby mode for up to 3 seconds.
- 2. The Shimmer unit will read the configuration file or infomem and create the required directories on the SD card.
- 3. The Shimmer unit will start logging.

Alternatively, if the Shimmer is powered on or reset while it **is** on the dock, the three steps above will begin as soon as the Shimmer is removed from the dock.

In either case, logging will continue until the Shimmer unit is reset, powered off, replaced in the dock or the battery runs out, whichever happens soonest.

Repeatedly resetting the Shimmer unit will result in multiple logging sessions on the SD card.

Warning: Please note that it is not recommend to dock the Shimmer unit while it is being configured. Ideally, the Shimmer unit should either be powered off or in standby mode whenever it is being placed on the dock.

User button (manual start)

This option relies on the orange user button on the *Shimmer3* enclosure to start/stop logging and is enabled if the "Push Button Start" option is enabled from ShimmerLog or ShimmerCapture (*i.e.*, configuration file line entry: "userbutton=1").

With this setting, undocking the Shimmer unit or powering it on off the dock will trigger the same steps 1 and 2 as above (*i.e.*, standby for 3 seconds followed by configuration). However, logging will not start until the user button is pressed by the user.



If the user button is pressed while the Shimmer unit is on the dock, it will have no effect on logging.

If the user button is pressed while the Shimmer unit is undocked, logging will start immediately.

Logging will continue until one of the following occurs: the user button is pressed, the Shimmer unit is reset, powered off, replaced in the dock, the battery runs out or the Shimmer is told to stop by a host application, whichever happens soonest.

Repeatedly pressing the user button and or docking the *Shimmer3* unit (to stop logging) will result in multiple logging sessions on the SD card.



3. Configuring the Shimmer3

3.1. Configuration and calibration loading priority

To use this firmware image, the user must provide the desired configuration parameters to each *Shimmer3*. For the *LogAndStream* firmware, configuration parameters are saved in two locations on the *Shimmer3*: a configuration file named *sdlog.cfg* located on the microSD card; and a flash memory location in the *Shimmer3*'s microcontroller called infomem. On initialisation, the *Shimmer3* will attempt to load configuration parameters from one of these sources in order of priority, initially from the configuration file and subsequently from infomem. If there is a problem loading parameters from either of the sources (*e.g.*, in the case of a new *Shimmer3*), the *Shimmer3* will revert to default values - as illustrated in Figure 3-1.




Figure 3-1 Shimmer3 configuration parameter and calibration value loading priority.

It is recommended to use the *ShimmerCapture or ShimmerLog* applications, available for download from the Shimmer website (www.shimmersensing.com), to avoid errors in writing the configuration file. However, it can be written in any text editor and saved to the SD card without the software application, by following the guidelines throughout this section.

The *LogAndStream* firmware is fully compatible with the *ShimmerLog* application. In this case, the user can configure the *Shimmer3* in the same manor as would previously have been done with the SDLog firmware, *i.e.*, with a *Shimmer Dock*. Please refer to the *ShimmerLog User Manual* for more information.



The ShimmerCapture application allows for configuration either through the Shimmer Dock or wirelessly over Bluetooth. In both cases, when the Shimmer3 is configured through the ShimmerCapture application, the new configuration parameters will be written to the Shimmer3's infomem. If the Shimmer3 is undocked when it is configured, the Shimmer3's microcontroller will automatically try to rewrite the configuration file with the new parameters. If the Shimmer3 is docked, as SD card access is handed over to the PC, the ShimmerCapture application will attempt to write the new configuration file.

Calibration parameters are not used in this firmware image to provide calibrated data, but they are saved in the configuration header of the raw data files on the microSD card. This is so that the data can calibrated post data collection. To supply calibration parameters, the user should create a folder called "Calibration" or "calibration" in the top-level SD card directory and should save the calibration file with the name "calibParams.ini" (name is case sensitive). The format of the file should match that output by the *Shimmer 9DOF Calibration Application* ².

If a calibration file is found by the firmware, the appropriate parameters from the file will be stored in the configuration header for the enabled sensors. As illustrated previously in Figure 3-1, if no calibration file is found, then the default calibration values will be stored for all sensors. If a calibration file is found but it does not contain parameters for a given enabled sensor, or for the configured sensor range, then default calibration values are stored for that sensor.

Note: When performing calibration with the *Shimmer 9DoF Calibration* application, the calibration data must be saved to the calibration file instead of using the "Save to Shimmer" option. The "Save to Shimmer" options only store the calibration data to the infomem of the *Shimmer3* and do not automatically update the calibration file. Therefore, as illustrated in Figure 3-1, when the *Shimmer3*'s power is cycled or the device is reset, the *LogAndStream* firmware will overwrite the 'new' calibration values in the infomem with the 'old' calibration file contents, if one is present.

3.2. Configuration file

The configuration file is read every time a new logging session is started on the Shimmer; thus, to change the configuration, just modify the *sdlog.cfg* file and reboot the Shimmer (there is no need to reprogram). The configuration file must be saved in the top level directory of the SD card (*i.e.*, not within a subfolder). This file allows configuration of all firmware features, as described throughout this section.

Each new line of the configuration file contains a single configuration parameter command. Users should note that the commands should be written **exactly** as they appear in the following sections so that they will be successfully parsed by the firmware. Any typographical errors, extra whitespaces, or repeated commands will cause parsing errors. It is the user's responsibility to ensure that the *sdlog.cfg* file is well-written. An example configuration file layout can be found in the appendices (*i.e.*, Section 6.2).

11

² For a more detailed description of IMU calibration parameters, refer to the *Shimmer 9DoF Calibration User Manual* and the *Shimmer IMU User Guide*, which are available for download from the members' area of www.shimmersensing.com.



3.3. Enabling sensors

The user may choose which sensors are enabled and disabled on the Shimmer for a particular experiment. The configuration file should contain a new line for each sensor enable/disable command. If any sensor is not included in the configuration file, it will be assumed to be disabled. The currently available sensors for the *Shimmer3* and their corresponding *sdlog.cfg* entries are listed in Table 3-1. There are currently no sensor conflicts for *Shimmer3*.

| Sensor | | Enable command | Disable command | |
|-----------------|-----------------|---|---|--|
| Low Noise Acce | lerometer | accel=1 | accel=0 | |
| Gyroscope | | gyro=1 | gyro=0 | |
| Magnetometer | | mag=1 | mag=0 | |
| Wide Range Acc | celerometer | accel_d=1 | accel_d=0 | |
| Battery voltage | | vbat=1 | vbat=0 | |
| External Expans | sion Channel 7 | extch7=1 | extch7=0 | |
| External Expans | sion Channel 6 | extch6=1 | extch6=0 | |
| External Expans | sion Channel 15 | extch5=1 | extch15=0 | |
| Internal Expans | ion Channel 1 | intch1=1 | intch1=0 | |
| Internal Expans | ion Channel 12 | intch12=1 | intch12=0 | |
| Internal Expans | ion Channel 13 | intch13=1 | intch13=0 | |
| Internal Expans | ion Channel 14 | intch14=1 | intch14=0 | |
| GSR | | gsr=1 | gsr=0 | |
| Pressure & Tem | perature | pres_bmp180=1 | pres_bmp180=0 | |
| Pulse/PPG | | intch13=1 exp power=1 | intch13=0 exp power=0 | |
| ECG | 24-bit | exg1_24bit=1 exg2_24bit=1 exp_power=1 | exg1_24bit=0 exg2_24bit=0 exp_Power=0 | |
| 16-bit | | exg1_16bit=1 exg2_16bit=1 exp_power=1 | exg1_16bit=0 exg2_16bit=0 exp_Power=0 | |
| EMG 24-bit | | exg1_24bit=1 | exg1_24bit=0 | |
| | 16-bit | exg1_16bit=1 | exg1_16bit=0 | |
| Strain Gauge | | str=1 | str=0 | |

Table 3-1 - Sensor enable/disable commands for Shimmer3. **Note:** the ECG and EMG channels can be set at either 16-bit or 24-bit but not both.

Note: To enable the ExG Test Signal to verify the correct operation of the *Shimmer3 ExG module* see the *ExG User Guide for ECG* or the *ExG User Guide for EMG,* both available for download on www.shimmersensing.com.



3.4. Configuring sensors

The following commands and options can be used to set sensor configuration parameters:

| Parameter | Example command | Valid options | Default |
|--------------------------|---------------------|----------------------------|----------------|
| Accel range | acc_range=0 | 0 (± 2.0 G) | 0 (± 2.0 G) |
| | | 1 (± 4.0 G) | |
| | | 2 (± 8.0 G) | |
| | | 3 (± 16.0 G) | |
| Gyro range | gyro_range=0 | 0 (±250 dps) | 0 (±250 dps) |
| | | 1 (±500 dps) | |
| | | 2 (±1000 dps) | |
| | | 3 (±2000 dps) | |
| Mag range | mg_range=1 | 1 (±1.3 Ga) | 1 (±1.3 Ga) |
| | | 2 (±1.9 Ga) | |
| | | 3 (±2.5 Ga) | |
| | | 4 (±4.0 Ga) | |
| | | 5 (±4.7 Ga) | |
| | | 6 (±5.6 Ga) | |
| | | 7 (±8.1 Ga) | |
| Accel internal data rate | acc_internal_rate=1 | 0 (power down) | 5 (100 Hz) |
| | | 1 (1.0 Hz) | |
| | | 2 (10 Hz) | |
| | | 3 (25 Hz) | |
| | | 4 (50 Hz) | |
| | | 5 (100 Hz) | |
| | | 6 (200 Hz) | |
| | | 7 (400 Hz) | |
| | | 8 (1.620KHz) | |
| | | 9 (1.344 KHz) | |
| Gyro internal data rate | gyro_samplingrate=0 | 0 – 255 (8000/value -1 Hz) | 155 (51.28 Hz) |
| Mag internal data rate | mg_internal_rate=0 | 0 (0.75 Hz) | 6 (75 Hz) |
| | | 1 (1.5 Hz) | |
| | | 2 (3.0 Hz) | |
| | | 3 (7.5 Hz) | |
| | | 4 (15.0 Hz) | |
| | | 5 (30.0 Hz) | |
| | | 6 (75.0 Hz) | |
| Accel Low Power Mode | acc_lpm=0 | 0 (disabled) | 0 (disabled) |
| | | 1 (enabled) | |
| Accel High Res Mode | acc_hrm=0 | 0 (disabled) | 0 (disabled) |
| | | 1 (enabled) | |
| GSR range | gs_range=0 | 0 (10kOhm – 56kOhm) | 4 (Auto Range) |
| | | 1 (56kOhm – 220kOhm) | |
| | | 2 (220kOhm – 680kOhm) | |
| | | 3 (680kOhm – 4.7MOhm) | |
| | | 4 (Auto Range) | |
| Pressure Resolution | pres_bmp180_prec=0 | 0 (Low) | 0 (Low) |
| | | 1 (Standard) | |
| | | 2 (High) | |
| | | 3 (Very High) | |

Table 3-2 - Sampling configuration options for Shimmer3.



Note: To change the gain and/or the data rate of the *ExG module*, see the *ExG User Guide for ECG* or the *ExG User Guide for EMG*, both available for download on www.shimmersensing.com.

3.5. Configuring an experiment

Logging preferences can be configured using the commands in Table 3-3. Further details describing these parameters are provided below the table.

| Parameter | Example command | Valid options | Default |
|------------------------|----------------------|--|---------|
| Sampling rate (in Hz) | sample_rate=51.2 | 0,, 1024 | 51.2 |
| Synchronisation | sync=0 | 0 (no time synchronisation) 1 (master-slave synchronisation) | 0 |
| Master | iammaster=0 | 0 (slave) 1 (master) | 0 |
| Broadcast interval (s) | interval=120 | integer value: 54,, 255 | 120 |
| User button | user_button_enable=0 | 0 (disabled) 1 (enabled) | 0 |
| Single-touch | singletouch=0 | 0 (disabled) 1 (enabled) | 0 |
| Number of Shimmers | Nshimmer=1 | integer value: 1,, 255 | 1 |
| Shimmer ID number | myid=1 | integer value: 1,, Nshimmer | 1 |
| Shimmer name | shimmername=shimmer1 | string with up to 11 characters | IDxxxx |
| Experiment ID | experimentid=expid | string with up to 11 characters | - |
| Configuration ID | configtime=1 | any 32-bit signed integer value | 0 |
| Master MAC Address | center=00066646b6af | MAC address of the Master BT Module | - |

Table 3-3 - Logging configuration options. **Note:** the multi-Shimmer sync features, as highlighted by the shaded rows, are not currently supported by the current version of LogAndStream.

Sampling rate

The Sampling rate is the frequency at which sampling should be performed. It should be noted that the actual sampling rate may not be exactly equal to this parameter because the firmware requires the sampling period to be an integer multiple of $1/32768 \, \text{s}$. For example, a desired sampling frequency of 500 Hz requires a sampling period of 2 ms which is equal to $65.536*(1/32768) \, \text{s}$. In the firmware, this sampling period will be rounded up to $66*(1/32768) = 2.014 \, \text{ms}$, so the true sampling rate will be $496.48 \, \text{Hz}$. The actual sampling rate that will be used (*i.e.*, the *True Fs*) can be calculated from the specified sampling rate (*Fs*) using the following formula:

$$True Fs = \frac{32768}{ceil \left(\frac{32768}{Desired Fs}\right)}$$

where the ceil() function means rounding up to the nearest integer.

Syncronisation, Master, Broadcast interval and Single-touch

As mentioned previously, the multi-Shimmer sync features (usually availble in the *Shimmer3 SDLog* and *BtStream* firmware), are not supported in the current version of *LogAndStream* firmware. For



compatibility reasons, the related configuration file entries (i.e., *Syncronisation, Master, Broadcast interval* and *Single-touch*) - as highlighted by the shaded rows in *Table 3-3*) are still present in the configuration file but currently have no effect.

User Button

The *User button* refers to the button on the baseboard which is accessible by pressing the circular orange button on the *Shimmer3* enclosure. If the button is enabled, then logging will begin when the user button is pressed. Logging will finish when **either** the user button is pressed again or the Shimmer is placed on the programming dock.

Note: the user button must not be enabled if the Heart Rate sensor is enabled due to conflicts caused by shared pins in hardware.

Number of Shimmers and Shimmer ID number

Number of Shimmers and *Shimmer ID number* are used to describe the total number of Shimmers in an experiment and the ID number of each individual Shimmer. Once again, these parameters do not have any explicit function in firmware.

Shimmer name

The *Shimmer name* parameter allows the user to specify a meaningful name for each device. It is used by the firmware to name the data directories on the SD card. For more information on the directory structure, see Section 5.

Experiment ID

The Experiment ID parameter allows the user to specify a meaningful name for an experiment. It is used by the firmware to name the data directories on the SD card. For more information on the directory structure, see Section 5.

Configuration ID

The Configuration ID is a numeric identifier for the experiment. It typically refers to the date and/or time at which the configuration file was created (but users may choose any valid numeric value that is meaningful to them). It is appended to the experiment ID name in the data directories to allow the user to identify different experiments in the case that the same experiment ID name is used multiple times. It is vital that this parameter is equal for all Shimmers in the experiment if synchronisation or single-touch start is enabled.

3.6. Infomem contents

The configuration parameter values are stored by the *Shimmer3* in the first 90 bytes of the infomem, which is the part of the *Shimmer3* memory that survives a reset or power cycle but is overwritten when the *Shimmer3* is reprogrammed. The format of the configuration data stored in infomem is as follows:



| Infomem Byte | Contents |
|--------------|--|
| 0 - 1 | Sampling rate |
| 2 | Buffer size |
| 3 - 5 | Selected sensors |
| 6 - 9 | Config bytes (Allows for 32 individual boolean settings) |
| 10 - 30 | Low Noise Accelerometer calibration values |
| 31 - 51 | Gyroscope calibration values |
| 52 - 72 | Magnetometer calibration values |
| 73 - 93 | Wide Range Accelerometer calibration values |

Table 3-4 Informem layout overview.

Selected Sensors - Infomem Bytes 3 to 5

The Selected sensors bytes have a single bit assigned to each sensor as follows:

| Bit | | Property | | | |
|----------------|---|------------------------------------|--|--|--|
| | 7 | Low Noise Accelerometer. | | | |
| m | 6 | Gyroscope. | | | |
| yte | 5 | Magnetometer. | | | |
| n B | 4 | ExG1_24BIT. | | | |
| ner | 3 | ExG2_24BIT. | | | |
| Infomem Byte 3 | 2 | GSR. | | | |
| 드 | 1 | External Expansion ADC Channel 7. | | | |
| | 0 | External Expansion ADC Channel 6. | | | |
| | 7 | Strain Gauge. | | | |
| 4 | 6 | Not yet assigned. | | | |
| yte | 5 | Battery Monitor. | | | |
| Infomem Byte 4 | 4 | Wide Range Accelerometer. | | | |
| ner | 3 | External Expansion ADC Channel 15. | | | |
| for | 2 | Internal Expansion ADC Channel 1. | | | |
| 드 | 1 | Internal Expansion ADC Channel 12. | | | |
| | 0 | Internal Expansion ADC Channel 13. | | | |
| | 7 | Internal Expansion ADC Channel 14. | | | |
| ιO | 6 | MPU9150 Accelerometer. | | | |
| yte | 5 | MPU9150 Magnetometer. | | | |
| n B | 4 | ExG1_16BIT. | | | |
| Infomem Byte 5 | 3 | ExG2_16BIT. | | | |
| for | 2 | BMP180 Pressure. | | | |
| 드 | 1 | BMP180 Temperature. | | | |
| | 0 | MSP430 Temperature. | | | |



Config Bytes - Infomem Bytes 6 to 9

The *Config bytes* contain the following parameters:

| Infomem B | Infomem Byte 6 - Config Setup Byte 0 | | | | |
|------------|---|--|--|--|--|
| Bits 7 – 4 | Wide Range (LSM303DLHC) Accelerometer Data Rate. | | | | |
| Bits 3 – 2 | Wide Range (LSM303DLHC) Accelerometer Range. | | | | |
| Bit 1 | Wide Range (LSM303DLHC) Accelerometer Low Power Mode. | | | | |
| Bit 0 | Wide Range (LSM303DLHC) Accelerometer High Resolution Mode. | | | | |
| Infomem B | yte 7 - Config Setup Byte 1 | | | | |
| Bits 7 – 0 | MPU9150 Data Rate. | | | | |
| Infomem B | yte 8 - Config Setup Byte 2 | | | | |
| Bits 7 – 5 | (LSM303DLHC) Magnetometer Range. | | | | |
| Bits 4 – 2 | (LSM303DLHC) Magnetometer Data Rate. | | | | |
| Bit 1 - 0 | MPU9150 Gyroscope Range. | | | | |
| Infomem B | yte 9 - Config Setup Byte 3 | | | | |
| Bits 7 – 6 | MPU9150 Accelerometer Range. | | | | |
| Bits 5 – 4 | BMP180 Pressure Resolution. | | | | |
| Bit 3 - 1 | GSR Range | | | | |
| Bit 0 | Internal Expansion Power Enable | | | | |

Calibration Parameters - Infomem Bytes 10 to 93

The calibration parameters for the inertial measurement units (accelerometer, gyroscope and magnetometer) consist of a three-element offset bias vector, a three-element sensitivity vector and a 3x3-element alignment matrix.³ The structure of these values when they are sent to/from the *Shimmer3* and stored in infomem is as follows:

- Each of the 3 offset bias vector values are stored as 16-bit signed integers (big endian) and are contained in bytes 0-5.
- Each of the 3 sensitivity vector values are stored as 16-bit signed integers (big endian) and are contained in bytes 6-11.
- Each of the 9 alignment matrix values are stored as 8-bit signed integers and are contained in bytes 12-20.

³ For a more detailed description of IMU calibration parameters, refer to the *Shimmer 9DoF Calibration User Manual* and the *Shimmer IMU User Guide*.



4. Bluetooth Streaming

The first byte of every packet received by the *Shimmer3* or the "host" is an identifier, telling the receiver what action to carry out or how to interpret the subsequent bytes. The full list of identifiers that are used to interface with the *BtStream* application, can be found in the header file, *Shimmer.h*, which can be found in the Appendix in Section 6.3 of this document (most recent version available online).

For every packet that the *Shimmer3* receives, it sends an acknowledgement message (ACK_COMMAND_PROCESSED) back to the host, to acknowledge receipt of the command.

4.1. Set Commands

The "SET" commands are used to set the values of all of the configurable parameters:

- Enabled sensors.
- Sampling rate.
- Accelerometer, gyroscope, magnetometer range.
- Accelerometer, gyroscope, magnetometer data rate.
- Battery monitoring.
- Calibration parameters for Accelerometers, Gyroscope, Magnetometer.
- Blink LED.

The packets sent between the Shimmer3 and the PC for a SET command are shown in Figure 4-1.




Figure 4-1 Packets sent for SET commands

These commands require that further data be received by the *Shimmer3* after the identifier byte. e.g. the SET_SAMPLING_RATE_COMMAND identifier must be followed by a one-byte value representing the sampling rate that the *Shimmer3* is to use. Another example is the SET_A_ACCEL_CALIBRATION_COMMAND identifier, which must be followed by 21 bytes representing the accelerometer calibration parameters.

4.2. Get Commands

The "GET" commands are requests for information and require that the *Shimmer3* sends data back to the host. The packets sent between the *Shimmer3* and the PC for a SET command are shown in Figure 4-2.

Figure 4-2 Packets sent for GET commands

On receipt of a GET command, the *Shimmer3* will send an acknowledgement message and, then, it will prepare and send a packet containing the appropriate response identifier byte, followed by the data that was requested.

For example, when the *Shimmer3* receives the GET_SAMPLING_RATE_COMMAND, it must send the current value of the sampling rate back to the host; the response packet will contain two bytes: the first byte will be the SAMPLING_RATE_REPONSE identifier and the second byte will be the sampling rate value.

Similarly, if the *Shimmer3* receives a GET_A_ACCEL_CALIBRATION_COMMAND, it will send a packet whose first byte is the A_ACCEL_CALIBRATION_RESPONSE identifier, followed by 21 bytes representing the accelerometer calibration parameters.

The INQUIRY_COMMAND

The INQUIRY_COMMAND is issued by the "host" when it wants to know the entire configuration of the *Shimmer3*, like what is the sampling rate, what is the buffer size, to which channel is each enabled sensor assigned, etc. In response to this command, the *Shimmer3* will send a packet back to the "host" with the structure shown in Table 4-1.

| Byte | 0 | 1-2 | 3-6 | 7 | 8 | 9 | 10 | ••• | х |
|-------|--------|----------|-----------|-------|--------|-------|-------|-----|-------|
| Value | Packet | Sampling | Config | Num | Buffer | Chan1 | Chan2 | | ChanX |
| | Туре | rate | Bytes 0-3 | Chans | size | _ = | | | |

Table 4-1 Inquiry response packet format

where the Packet Type = INQUIRY_RESPONSE and the value in the channel fields (Chan1, Chan2, ..., ChanX) indicate exactly what data from which sensor will be contained in the equivalent field of the data packet. The total number of bytes sent by the *Shimmer3* will depend on how many data channels are active (i.e. which sensors are enabled).

Signal name, byte values and datatypes

Table 4-2 lists the values in the channel contents bytes of the Inquiry response packet along with the signal names and datatypes for the equivalent sensor signals (* in the datatype column denotes MSB first; otherwise LSB first).

| Signal Name | Byte Value | Signal Datatype |
|----------------------------|------------|-----------------|
| Low Noise Accelerometer X* | 0 | u12 |
| Low Noise Accelerometer Y* | 1 | u12 |
| Low Noise Accelerometer Z* | 2 | u12 |



| Battery | 3 | u12 |
|-----------------------------|----|------|
| Wide Noise Accelerometer X* | 4 | i16 |
| Wide Noise Accelerometer Y* | 5 | i16 |
| Wide Noise Accelerometer Z* | 6 | i16 |
| Gyroscope X* | 7 | i16* |
| Gyroscope Y* | 8 | i16* |
| Gyroscope Z* | 9 | i16* |
| Magnetometer X* | Α | i16* |
| Magnetometer Y* | В | i16* |
| Magnetometer Z* | С | i16* |
| External ADC 7 | D | u12 |
| External ADC 6 | E | u12 |
| External ADC 15 | F | u12 |
| Internal ADC 1 | 10 | u12 |
| Internal ADC 12 | 11 | u12 |
| Internal ADC 13 | 12 | u12 |
| Internal ADC 14 | 13 | u12 |
| BMP180 Temperature* | 1A | u16* |
| BMP180 Pressure* | 1B | u24* |
| GSR Raw | 1C | u16 |
| ExG_ADS1292R_1_STATUS | 1D | u8 |
| ExG_ADS1292R_1_CH1_24BIT | 1E | i24* |
| ExG_ADS1292R_1_CH2_24BIT | 1F | i24* |
| ExG_ADS1292R_2_STATUS | 20 | u8 |
| ExG_ADS1292R_2_CH1_24BIT | 21 | i24* |
| ExG_ADS1292R_2_CH2_24BIT | 22 | i24* |
| ExG_ADS1292R_1_CH1_16BIT | 23 | i16* |
| ExG_ADS1292R_1_CH2_16BIT | 24 | i16* |
| ExG_ADS1292R_2_CH1_16BIT | 25 | i16* |
| ExG_ADS1292R_2_CH2_16BIT | 26 | i16* |
| Strain Gauge High | 27 | u12 |
| Strain Gauge Low | 28 | u12 |

Table 4-2 Signal names, channel contents byte values and datatypes for available sensor signals

4.3. Action Commands

There are a number of available "ACTION" commands, which do not require that parameter values be sent between the PC and the *Shimmer3* but, instead, tell the *Shimmer3* what action it is to carry out. These include the START_STREAMING_COMMAND and STOP_STREAMING_COMMAND and the TOGGLE_LED_COMMAND.

4.4. Streaming

When the START_STREAMING_COMMAND is received by the *Shimmer3*, it will send an acknowledge message back to the "host" and start sampling sensor data. As the sensor data is sampled, the *Shimmer3* will prepare data packets and send them to the "host" over Bluetooth.

The *Buffer size* parameter determines the number of samples that are sent together in a single data packet. The structure of the data packet with Buffer size = 2 is shown in Table 4-3, where Packet Type = DATA_PACKET, TS denotes "Timestamp" and Ch denotes "Channel".



| Byte | 0 | 1 - 2 | 3 - 4 | 5 - 6 | (x-1) - x | (x+1) - (x+2) | (x+3) - (x+4) | (x+5) - (x+6) | | (2x-1) - 2x |
|-------|----------|-------|-------|-------|---------------|---------------|---------------|---------------|--|-------------|
| | Packet | | | | | | | | | |
| Value | Туре | TS | Ch1 | Ch2 | ChX | TS | Ch1 | Ch2 | | ChX |
| | Sample 1 | | | | Sample 2 | | | | | |

Table 4-3 Data packet structure (Buffer size =2)

If *Buffer size* were equal to 1, then the data packet would contain only one timestamp and one sample from each channel (i.e. the bytes denoted "Sample 1" in Table 4-3. If *Buffer size* were any integer value greater than 2, then subsequent timestamps and sample values for each channel would be appended at the end of the packet until the number of samples equals the buffer size.

Sensor data will continue to be sampled and streamed until a STOP_STREAMING_COMMAND is received by the *Shimmer3*.

By default, the application will sample the 3-axis accelerometer at 51.2 Hz and send the data using a data buffer of size 1.

4.5. Bluetooth Latency

Our lab tests have shown up to 100 ms of latency with considerable variation (> 50 ms). These measures result from multiple FIFOs in the data path, as expected in wireless data acquisition systems using conventional computing devices for the data end-points. Actual performance is strongly impacted by end-point system configuration and load.



5. Reading the SD card Data

It is recommended that the *ShimmerLog* software, which can be downloaded from the Shimmer website (www.shimmersensing.com) be used to read the data from the SD card. However, using the following guidelines, the user may use their platform of choice to read and parse the data.

5.1. SD card directory structure

The data directory structure on the SD card is as follows:

- A folder called "data" is created by the firmware in the SD card top-level directory if it does not already exist.
- A subfolder is created within the *data* folder each time a new "experiment id" is encountered in the *sdlog.cfq* file (default is "default exp").
- A new subfolder is created within the appropriate experiment folder each time logging starts. The naming convention for this folder is the Shimmer name specified in the *sdlog.cfg* file (default is the abbreviated form of the *Shimmer3* id number (e.g. *IDbf5e*)), followed by a three digit number which is sequentially incremented at new logging session.
- For example, the first time a *Shimmer3* with name "device1" in experiment "experiment1" logs data to the SD card it will create the folder data/experiment1/device1-000/. The next time it stops and subsequently restarts logging, it will create data/experiment1/device1-001/, etc.
- Within this subfolder, data is logged in files which are named sequentially with a three digit number indicating the order of logging, starting with 000. The file is closed after 1 hour of continuous data logging and a new file is opened in the same subfolder, such that the second file is called 001 and so on. The last file is closed when logging stops and a new file (in a different folder) is created the next time logging begins.

5.2. Parsing the raw data file

Raw data file configuration header

Each RAW data file (e.g. data/experiment1/device1-000/000), begins with a header which contains the Shimmer3 configuration information. The structure of this header is outlined in the tables below, where each row of each table represents one byte and the individual bits are separated, where appropriate, depending on whether each bit represents an independent binary value or the entire byte contains a single value. Empty cells and cells with constant values represent bits that are reserved for future use. For Shimmer3, there are a total 256 bytes in the header.

Byte # 0 - 9: Enabled Sensors

| Byte | Bit 7 | Bit 6 | Bit 5 | Bit 4 | Bit 3 | Bit 2 | Bit 1 | Bit 0 | |
|------|-----------------------|-----------------------|-------|-------|-------|-------|-------|-------|--|
| 0 | ADC Sample Rate (LSB) | | | | | | | | |
| 1 | | ADC Sample Rate (MSB) | | | | | | | |
| 2 | | | | | | | | | |



| 3 | Acc (LN)* | Gyr | Mag | EXG1_24BI T | EXG2_24BI T | GSR | ExtCh 7 | ExtCh 6 |
|---|-----------|-------------------|-----------------|----------------|----------------|---------|-----------|-----------|
| 4 | strain | | Battery | Acc (WR)* | ExtCh 15 | IntCh 1 | IntCh 12 | IntCh13 |
| 5 | IntCh 14 | MPU9150_ ACCEL | MPU9150_ MAG | EXG1_16BI T | EXG2_16BI T | pres | | |
| 6 | | | | | | | | |
| 7 | | | | | | | | |
| 8 | | LSM303 Digt | ial Accel Rate | | D Acce | l Range | Accel LPM | Accel HRM |
| 9 | | | | MPU9150 | Gyro Rate | | • | • |

^{*} There are multiple accelerometers on the *Shimmer3* base board. The low noise (LN) analog accelerometer on the KXRB5-2042 chip is enabled via Byte 3, Bit 7 (Acc (LN)), whilst the wide range digital accelerometer on the LSM303DLHC chip is enabled via Byte 4, Bit 4 (Acc (WR)).

Byte # 10 - 19: Trial Configuration

| Byte | Bit 7 | Bit 6 | Bit 5 | Bit 4 | Bit 3 | Bit 2 | Bit 1 | Bit 0 |
|------|---------|--------------|--------------------|--------------------|--------------|--------------------|--------|---------|
| 10 | LSI | M303 Mag Ran | ige | LS | SM303 Mag Ra | MPU9150 Gyro Range | | |
| 11 | MPU_ACC | EL_RANGE | Pressure Precision | | | GSR Range | | EXP_PWR |
| 12 | | | | | | | | |
| 13 | | | | | | | | |
| 14 | | | | | | | | |
| 15 | | | | | | | | |
| 16 | | | User | | 1 | Sync | Master | 0 |
| | | | button | | | | | |
| 17 | Single | Accel LPM | Accel HRM | txco | | | | |
| | touch | | | | | | | |
| 18 | | | | Broadcast interval | | | | |
| 19 | | | | | | | | |

Byte # 30 – 39: Firmware and Shimmer Parameters

| Byte | Bit 7 | Bit 6 | Bit 5 | Bit 4 | Bit 3 | Bit 2 | Bit 1 | Bit 0 | | | | | |
|------|-------|----------------------|-------|------------|-----------------|-------|-------|-------|--|--|--|--|--|
| 30 | | | | Shimmer | Version (MSB) | ' | | • | | | | | |
| 31 | | ShimmerVersion (LSB) | | | | | | | | | | | |
| 32 | | MyTriald | | | | | | | | | | | |
| 33 | | Nshimmer | | | | | | | | | | | |
| 34 | | | | FW Version | on - Type (MSB) | | | | | | | | |
| 35 | | | | FW Versi | on - Type (LSB) | | | | | | | | |
| 36 | | | | FW Versio | n - Major (MSB | 3) | | | | | | | |
| 37 | | | | FW Versi | on - Major (LSB | | | | | | | | |
| 38 | | | | FW Ve | rsion - Minor | | | | | | | | |
| 39 | | | | FW Ver | sion - Release | | | | | | | | |

Shimmer Version indicates the hardware version for which the code was compiled. Shimmer3 is denoted by 3.

The FW Version bytes define the firmware version. The Type field will always be 2 for SDLog firmware images. Major and Minor versions are indicated by a two- and one-byte value, respectively



(e.g. for LogAndStreamg v1.0, Major = 1 and Minor = 0). The Release field can be ignored by users; it will have a value of 0.

Byte # 40 - 51: Reserved for future use

Byte # 52 - 55: Configuration Time

| Byte | Bit 7 | Bit 6 | Bit 5 | Bit 4 | Bit 3 | Bit 2 | Bit 1 | Bit 0 | | | |
|------|---------------|---------------------|-------|------------|------------|-------|-------|-------|--|--|--|
| 52 | | Config Time 0 (MSB) | | | | | | | | | |
| 53 | Config Time 1 | | | | | | | | | | |
| 54 | | Config Time 2 | | | | | | | | | |
| 55 | | | | Config Tir | ne 3 (LSB) | | | | | | |

This value comes from the *Configuration ID* parameter from Section 3.5.

Byte # 56-181: Calibration parameters

Bytes 56 to 181 store the calibration data for individual *Shimmer3* sensors. These parameters are used to calibrate the RAW data during post-experimentation analysis. The layout and description of these bytes can be found in the Appendix section of this document (*i.e.*, Section 6.2).

Byte # 182 - 251: Reserved for future use

Byte # 252 - 255: Initial Timestamp

| Byte | Bit 7 | Bit 6 | Bit 5 | Bit 4 | Bit 3 | Bit 2 | Bit 1 | Bit 0 | | |
|------|--------------------|--------------------|-------|-------------|----------|-------|-------|-------|--|--|
| 252 | Initial time stamp | | | | | | | | | |
| 253 | Initial time stamp | | | | | | | | | |
| 254 | | Initial time stamp | | | | | | | | |
| 255 | | | | Initial tin | ne stamp | | | | | |

Immediately after the configuration header, there is a 4-byte timestamp, whose value is the time on the local *Shimmer3* clock (in units of ticks of the 32 kHz clock) when the first sample in that file was recorded. It is used for aligning the data from multiple Shimmer devices to a common clock in multi-shimmer sync applications. As the multi-shimmer sync feature is not available in the current LogAndStream firmware, these bytes can be ignored.

Raw data file format

Following these 256 bytes, the logged sensor data begins (i.e. at the 257th byte in the file). The data is written in blocks of up to 512 bytes. These blocks are continuously written until data has been logged to a given file for one hour, at which time the current file is closed and a new file is opened, with the same format as the current one (i.e. the configuration header is repeated in every file).

The exact number of bytes per block and the interpretation of those bytes of data depend on the enabled sensors. The total number of bytes written per block, *Bp*, is given by:

$$Bp = N \times (Bs + 2)$$



where:

Bs = Nc3 * 3 + Nc2 * 2 + Nc1, is the number of bytes per sample, which depends on the number of 3byte channels (Nc3), number of 2byte channels (Nc2) and the number of 1-byte channels (Nc1) required by the enabled sensors;

N = floor(512/(Bs + 2)), is the integer number of samples per block, with a buffer of 512 bytes and the extra 2 bytes per sample are the 16-bit timestamps from the *Shimmer3*'s local 32 kHz clock.

All of the bytes in the block are interleaved local 16-bit timestamps and the associated sensor samples. For example, with the accelerometer and the gyroscope enabled, there are three 2-byte channels for the accelerometer and three 2-byte channels for the gyroscope in each sample and the data would be as shown in the following table until a total of Bp = 504 bytes is reached (Bs = 12; N = $\frac{1}{1000}$) floor((512)/(12 + 2)) = 36).

| Byte | Parameter | Sensor |
|------|---------------------------------|----------|
| 0 | timestamp 0 byte 0 (MSB) | |
| 1 | timestamp 0 byte 1 (LSB) | |
| 2 | sample 0 channel 0 byte 0 (MSB) | (XAccel) |
| 3 | sample 0 channel 0 byte 1 (LSB) | (AACCEI) |
| 4 | sample 0 channel 1 byte 0 (MSB) | (YAccel) |
| 5 | sample 0 channel 1 byte 1 (LSB) | (TACCEI) |
| 6 | sample 0 channel 2 byte 0 (MSB) | (ZAccel) |
| 7 | sample 0 channel 2 byte 1 (LSB) | (ZACCEI) |
| 8 | sample 0 channel 3 byte 0 (MSB) | (XGyro) |
| 9 | sample 0 channel 3 byte 1 (LSB) | (Adylo) |
| 10 | sample 0 channel 4 byte 0 (MSB) | (YGyro) |
| 11 | sample 0 channel 4 byte 1 (LSB) | (Taylo) |
| 12 | sample 0 channel 5 byte 0 (MSB) | (ZGyro) |
| 13 | sample 0 channel 5 byte 1 (LSB) | (ZGylo) |
| 14 | timestamp 1 byte 0 (MSB) | |
| 15 | timestamp 1 byte 1 (LSB) | |
| 16 | sample 1 channel 0 byte 0 (MSB) | (XAccel) |
| 17 | sample 1 channel 0 byte 1 (LSB) | (AACCEI) |
| | | |

Table 5-1 - Example of the raw data format for a file stored to the SD card of a Shimmer3.

Order of raw data channels

The order in which the sensors appear in the data channels depends on which sensors are enabled. The firmware will assign the channels in the order outlined in Table 5-2 (top to bottom). Note that "Accel (LN)" refers to the low noise accelerometer, whilst "Accel (WR)" refers to the wide range accelerometer. The number of channels per sensor is listed beside the channel name. The endianess of the bytes for each channel is also specified.

| Channel Channel cont | nts Number of | Bytes per | Endian |
|----------------------|---------------|-----------|--------|
|----------------------|---------------|-----------|--------|



| type | | channels | channel | |
|----------|---------------------|----------|---------|--------|
| Analog | Accel (LN) | 3 | 2 | little |
| channels | Battery | 1 | 2 | little |
| | Ext Exp A7 | 1 | 2 | little |
| | Ext Exp A6 | 1 | 2 | little |
| | Ext Exp A15 | 1 | 2 | little |
| | Int Exp A12 | 1 | 2 | little |
| | Int Exp A13 | 1 | 2 | little |
| | Int Exp A14 | 1 | 2 | little |
| | Strain High | 1 | 2 | little |
| | Strain Low | 1 | 2 | little |
| | Int Exp A1 | 1 | 2 | little |
| Digital | Gyro_mpu | 3 | 2 | big |
| channels | Accel_lsm (WR) | 3 | 2 | little |
| | Mag_lsm | 3 | 2 | big |
| | Accel_mpu | 3 | 2 | big |
| | Mag_mpu | 3 | 2 | little |
| | Temperature_bmp180 | 1 | 2 | big |
| | Pressure_bmp180 | 1 | 3 | big |
| | EXG1(24-bit/16-bit) | 1 | 7/5 | big |
| | EXG2(24-bit/16-bit) | 1 | 7/5 | big |

Table 5-2 - Order of data channels for Shimmer3 (top to bottom).



6. Appendices

6.1. Example of the Shimmer3 configuration file (sdlog.cfg)

```
accel=1
gyro=1
mag=1
accel d=1
vbat=1
extch7=0
extch6=0
extch15=0
intch1=0
intch12=0
intch13=0
intch14=0
gsr=0
str=0
pres bmp180=0
exg1 24bit=0
exg2_24bit=0
exg1_16bit=0
exg2_16bit=0
EXG ADS1292R 1 CONFIG1=0
EXG ADS1292R 1 CONFIG2=163
EXG ADS1292R 1 LOFF=16
EXG ADS1292R 1 CH1SET=69
EXG ADS1292R 1 CH2SET=5
EXG ADS1292R 1 RLD SENS=0
EXG ADS1292R 1 LOFF SENS=0
EXG ADS1292R 1 LOFF STAT=0
EXG ADS1292R 1 RESP1=2
EXG ADS1292R 1 RESP2=1
EXG ADS1292R 2 CONFIG1=2
EXG ADS1292R 2 CONFIG2=163
EXG ADS1292R 2 LOFF=16
EXG ADS1292R 2 CH1SET=5
EXG_ADS1292R 2 CH2SET=5
EXG ADS1292R 2 RLD SENS=0
EXG_ADS1292R_2_LOFF_SENS=0
EXG ADS1292R 2 LOFF STAT=0
EXG ADS1292R 2 RESP1=2
EXG ADS1292R 2 RESP2=1
exp_power=0
acc range=1
acc internal rate=5
acc_lpm=0
acc hrm=0
gyro samplingrate=155
gyro range=1
mg internal rate=6
mg range=1
user button enable=1
```



sample_rate=51.2
iammaster=0
sync=1
interval=120
singletouch=0
center=00066646b6af
myid=1
Nshimmer=2
shimmername=device1
experimentid=expidname
configtime=1234567

6.2. SD data file configuration header - sensor calibration bytes

Byte # 56-76: ExG Calibration

For *Shimmer3*, these parameters refer to the ExG module which is made up of two ADS1292R chips.

| 56 | | N 13 | | | | | | | | | |
|----|---------------------------|---------------------------|--------|--------|--------|--------|-----|--|--|--|--|
| | | NV_EXG_ADS1292R_1_CONFIG1 | | | | | | | | | |
| 57 | NV_EXG_ADS1292R_1_CONFIG2 | | | | | | | | | | |
| 58 | | 1 | NV_EXC | G_ADS1 | .292R_ | 1_LOFF | : | | | | |
| 59 | | N | V_EXG | _ADS12 | 92R_1 | _CH1SI | ΞT | | | | |
| 60 | | N | V_EXG | _ADS12 | 92R_1 | _CH2SI | ΞT | | | | |
| 61 | | NV | _EXG_A | ADS129 | 2R_1_ | RLD_SE | NS | | | | |
| 62 | | NV_ | EXG_A | DS129 | 2R_1_l | .OFF_S | ENS | | | | |
| 63 | | NV_ | EXG_A | DS129 | 2R_1_l | .OFF_S | TAT | | | | |
| 64 | | N | V_EXG | _ADS1 | 292R_1 | L_RESP | 1 | | | | |
| 65 | NV_EXG_ADS1292R_1_RESP2 | | | | | | | | | | |
| 66 | | NV | _EXG_ | ADS12 | 92R_2_ | CONFI | G1 | | | | |
| 67 | п | NV | _EXG_ | ADS12 | 92R_2_ | CONFI | G2 | | | | |
| 68 | | ſ | NV_EXC | G_ADS1 | .292R_ | 2_LOFF | | | | | |
| 69 | | N | V_EXG_ | ADS12 | 92R_2 | _CH1SI | ΞT | | | | |
| 70 | | N | V_EXG | _ADS12 | 92R_2 | _CH2SI | ΞT | | | | |
| 71 | | NV. | _EXG_/ | ADS129 | 2R_2_ | RLD_SE | NS | | | | |
| 72 | | NV_ | EXG_A | DS129 | 2R_2_L | .OFF_S | ENS | | | | |
| 73 | | NV_ | EXG_A | DS129 | 2R_2_l | .OFF_S | TAT | | | | |
| 74 | | N | V_EXG | _ADS1 | 292R_2 | 2_RESP | 1 | | | | |
| 75 | | N | V_EXG | _ADS1 | 292R_2 | 2_RESP | 2 | | | | |

Byte # 76-96: Digital Accelerometer Calibration

For *Shimmer3*, these parameters refer to the LSM303DLHC accelerometer. There are other accelerometers also populated on the *Shimmer3* circuit board; the associated calibration parameters for the other devices are stored in a later section of the configuration header.

| Byte | Bit 7 | Bit 6 | Bit 5 | Bit 4 | Bit 3 | Bit 2 | Bit 1 | Bit 0 | | |
|------|-------|--|----------|---------|---------|---------|-------|-------|--|--|
| 76 | | Digital Accel Calibration Offset X (MSB) | | | | | | | | |
| 77 | | Digita | al Accel | Calibra | ation O | ffset X | (LSB) | | | |



| 78Digital Accel Calibration Offset Y (MSB)79Digital Accel Calibration Offset Y (LSB)80Digital Accel Calibration Offset Z (MSB)81Digital Accel Calibration Offset Z (LSB)82Digital Accel Calibration Gain X (MSB)83Digital Accel Calibration Gain Y (MSB)84Digital Accel Calibration Gain Y (MSB)85Digital Accel Calibration Gain Z (MSB)86Digital Accel Calibration Gain Z (MSB)87Digital Accel Calibration Gain Z (LSB)88Digital Accel Calibration Align XX89Digital Accel Calibration Align XZ90Digital Accel Calibration Align XZ91Digital Accel Calibration Align YX92Digital Accel Calibration Align YY93Digital Accel Calibration Align YZ94Digital Accel Calibration Align ZX95Digital Accel Calibration Align ZY | | |
|---|----|--|
| B0 Digital Accel Calibration Offset Z (MSB) B1 Digital Accel Calibration Offset Z (LSB) B2 Digital Accel Calibration Gain X (MSB) B3 Digital Accel Calibration Gain X (LSB) B4 Digital Accel Calibration Gain Y (MSB) B5 Digital Accel Calibration Gain Y (LSB) B6 Digital Accel Calibration Gain Z (MSB) B7 Digital Accel Calibration Gain Z (MSB) B8 Digital Accel Calibration Align XX B9 Digital Accel Calibration Align XX 90 Digital Accel Calibration Align XZ 91 Digital Accel Calibration Align YX 92 Digital Accel Calibration Align YY 93 Digital Accel Calibration Align YZ 94 Digital Accel Calibration Align YZ 95 Digital Accel Calibration Align YZ 96 Digital Accel Calibration Align YZ 97 Digital Accel Calibration Align YZ 98 Digital Accel Calibration Align YZ | 78 | Digital Accel Calibration Offset Y (MSB) |
| B1 Digital Accel Calibration Offset Z (LSB) B2 Digital Accel Calibration Gain X (MSB) B3 Digital Accel Calibration Gain X (LSB) B4 Digital Accel Calibration Gain Y (MSB) B5 Digital Accel Calibration Gain Y (LSB) B6 Digital Accel Calibration Gain Z (MSB) B7 Digital Accel Calibration Gain Z (LSB) B8 Digital Accel Calibration Align XX B9 Digital Accel Calibration Align XY 90 Digital Accel Calibration Align XZ 91 Digital Accel Calibration Align YX 92 Digital Accel Calibration Align YY 93 Digital Accel Calibration Align YZ 94 Digital Accel Calibration Align YZ | 79 | Digital Accel Calibration Offset Y (LSB) |
| B2 Digital Accel Calibration Gain X (MSB) B3 Digital Accel Calibration Gain X (LSB) B4 Digital Accel Calibration Gain Y (MSB) B5 Digital Accel Calibration Gain Y (LSB) B6 Digital Accel Calibration Gain Z (MSB) B7 Digital Accel Calibration Gain Z (LSB) B8 Digital Accel Calibration Align XX B9 Digital Accel Calibration Align XY 90 Digital Accel Calibration Align XZ 91 Digital Accel Calibration Align YX 92 Digital Accel Calibration Align YY 93 Digital Accel Calibration Align YZ 94 Digital Accel Calibration Align YZ | 80 | Digital Accel Calibration Offset Z (MSB) |
| B3 Digital Accel Calibration Gain X (LSB) B4 Digital Accel Calibration Gain Y (MSB) B5 Digital Accel Calibration Gain Y (LSB) B6 Digital Accel Calibration Gain Z (MSB) B7 Digital Accel Calibration Gain Z (LSB) B8 Digital Accel Calibration Align XX B9 Digital Accel Calibration Align XY 90 Digital Accel Calibration Align XZ 91 Digital Accel Calibration Align YX 92 Digital Accel Calibration Align YY 93 Digital Accel Calibration Align YZ 94 Digital Accel Calibration Align ZX | 81 | Digital Accel Calibration Offset Z (LSB) |
| B4 Digital Accel Calibration Gain Y (MSB) B5 Digital Accel Calibration Gain Y (LSB) B6 Digital Accel Calibration Gain Z (MSB) B7 Digital Accel Calibration Gain Z (LSB) B8 Digital Accel Calibration Align XX B9 Digital Accel Calibration Align XY 90 Digital Accel Calibration Align XZ 91 Digital Accel Calibration Align YX 92 Digital Accel Calibration Align YY 93 Digital Accel Calibration Align YZ 94 Digital Accel Calibration Align ZX | 82 | Digital Accel Calibration Gain X (MSB) |
| B5 Digital Accel Calibration Gain Y (LSB) B6 Digital Accel Calibration Gain Z (MSB) B7 Digital Accel Calibration Gain Z (LSB) B8 Digital Accel Calibration Align XX B9 Digital Accel Calibration Align XY 90 Digital Accel Calibration Align XZ 91 Digital Accel Calibration Align YX 92 Digital Accel Calibration Align YY 93 Digital Accel Calibration Align YZ 94 Digital Accel Calibration Align ZX | 83 | Digital Accel Calibration Gain X (LSB) |
| B6 Digital Accel Calibration Gain Z (MSB) B7 Digital Accel Calibration Gain Z (LSB) B8 Digital Accel Calibration Align XX B9 Digital Accel Calibration Align XY 90 Digital Accel Calibration Align XZ 91 Digital Accel Calibration Align YX 92 Digital Accel Calibration Align YY 93 Digital Accel Calibration Align YZ 94 Digital Accel Calibration Align ZX | 84 | Digital Accel Calibration Gain Y (MSB) |
| 87 Digital Accel Calibration Gain Z (LSB) 88 Digital Accel Calibration Align XX 89 Digital Accel Calibration Align XY 90 Digital Accel Calibration Align XZ 91 Digital Accel Calibration Align YX 92 Digital Accel Calibration Align YY 93 Digital Accel Calibration Align YZ 94 Digital Accel Calibration Align ZX | 85 | Digital Accel Calibration Gain Y (LSB) |
| 88 Digital Accel Calibration Align XX 89 Digital Accel Calibration Align XY 90 Digital Accel Calibration Align XZ 91 Digital Accel Calibration Align YX 92 Digital Accel Calibration Align YY 93 Digital Accel Calibration Align YZ 94 Digital Accel Calibration Align ZX | 86 | Digital Accel Calibration Gain Z (MSB) |
| 89 Digital Accel Calibration Align XY 90 Digital Accel Calibration Align XZ 91 Digital Accel Calibration Align YX 92 Digital Accel Calibration Align YY 93 Digital Accel Calibration Align YZ 94 Digital Accel Calibration Align ZX | 87 | Digital Accel Calibration Gain Z (LSB) |
| 90 Digital Accel Calibration Align XZ 91 Digital Accel Calibration Align YX 92 Digital Accel Calibration Align YY 93 Digital Accel Calibration Align YZ 94 Digital Accel Calibration Align ZX | 88 | Digital Accel Calibration Align XX |
| 91 Digital Accel Calibration Align YX 92 Digital Accel Calibration Align YY 93 Digital Accel Calibration Align YZ 94 Digital Accel Calibration Align ZX | 89 | Digital Accel Calibration Align XY |
| 92 Digital Accel Calibration Align YY 93 Digital Accel Calibration Align YZ 94 Digital Accel Calibration Align ZX | 90 | Digital Accel Calibration Align XZ |
| 93 Digital Accel Calibration Align YZ 94 Digital Accel Calibration Align ZX | 91 | Digital Accel Calibration Align YX |
| 94 Digital Accel Calibration Align ZX | 92 | Digital Accel Calibration Align YY |
| 0 0 | 93 | Digital Accel Calibration Align YZ |
| 95 Digital Accel Calibration Align ZY | 94 | Digital Accel Calibration Align ZX |
| | 95 | Digital Accel Calibration Align ZY |
| 96 Digital Accel Calibration Align ZZ | 96 | Digital Accel Calibration Align ZZ |

Byte # 97 - 117: Gyroscope Calibration

| Byte | Bit 7 | Bit 6 | Bit 5 | Bit 4 | Bit 3 | Bit 2 | Bit 1 | Bit 0 | | | |
|------|---------------------------------|---------------------------------|----------|----------|----------|----------|-------|-------|--|--|--|
| 97 | Gyro Calibration Offset X (MSB) | | | | | | | | | | |
| 98 | Gyro Calibration Offset X (LSB) | | | | | | | | | | |
| 99 | Gyro Calibration Offset Y (MSB) | | | | | | | | | | |
| 100 | | Gyro Calibration Offset Y (LSB) | | | | | | | | | |
| 101 | | Gy | yro Cali | bration | offset | t Z (MS | В) | | | | |
| 102 | | G | yro Cal | ibratio | n Offse | t Z (LSE | 3) | | | | |
| 103 | | G | iyro Ca | libratio | n Gain | X (MSE | 3) | | | | |
| 104 | | (| Gyro Ca | libratio | n Gain | X (LSB |) | | | | |
| 105 | | G | iyro Ca | libratio | n Gain | Y (MSE | 3) | | | | |
| 106 | | (| Gyro Ca | libratio | n Gain | Y (LSB |) | | | | |
| 107 | | G | iyro Ca | libratio | n Gain | Z (MSE | 3) | | | | |
| 108 | | (| Gyro Ca | libratio | n Gain | Z (LSB |) | | | | |
| 109 | | | Gyro | Calibra | tion Ali | gn XX | | | | | |
| 110 | | | Gyro | Calibra | tion Ali | gn XY | | | | | |
| 111 | | | Gyro | Calibra | tion Ali | gn XZ | | | | | |
| 112 | Gyro Calibration Align YX | | | | | | | | | | |
| 113 | Gyro Calibration Align YY | | | | | | | | | | |
| 114 | Gyro Calibration Align YZ | | | | | | | | | | |
| 115 | Gyro Calibration Align ZX | | | | | | | | | | |
| 116 | | Gyro Calibration Align ZY | | | | | | | | | |
| 117 | | | Gyro | Calibra | tion Ali | gn ZZ | | | | | |

Byte # 118 - 138: Magnetometer calibration

| Byte | Bit 7 | Bit 6 | Bit 5 | Bit 4 | Bit 3 | Bit 2 | Bit 1 | Bit 0 |
|------|-------|-------|-------|-------|-------|-------|-------|-------|



| 118 | Mag Calibration Offset X (MSB) |
|-----|--------------------------------|
| 119 | Mag Calibration Offset X (LSB) |
| 120 | Mag Calibration Offset Y (MSB) |
| 121 | Mag Calibration Offset Y (LSB) |
| 122 | Mag Calibration Offset Z (MSB) |
| 123 | Mag Calibration Offset Z (LSB) |
| 124 | Mag Calibration Gain X (MSB) |
| 125 | Mag Calibration Gain X (LSB) |
| 126 | Mag Calibration Gain Y (MSB) |
| 127 | Mag Calibration Gain Y (LSB) |
| 128 | Mag Calibration Gain Z (MSB) |
| 129 | Mag Calibration Gain Z (LSB) |
| 130 | Mag Calibration Align XX |
| 131 | Mag Calibration Align XY |
| 132 | Mag Calibration Align XZ |
| 133 | Mag Calibration Align YX |
| 134 | Mag Calibration Align YY |
| 135 | Mag Calibration Align YZ |
| 136 | Mag Calibration Align ZX |
| 137 | Mag Calibration Align ZY |
| 138 | Mag Calibration Align ZZ |

Byte # 139 - 159: Analog Accelerometer Calibration

| Byte | Bit 7 Bit 6 Bit 5 Bit 4 Bit 3 Bit 2 Bit 1 Bit 0 | | | | | | | | | |
|------|---|---|---------|----------|----------|----------|-------|--|--|--|
| 139 | Analog Accel Calibration Offset X (MSB) | | | | | | | | | |
| 140 | Analog Accel Calibration Offset X (LSB) | | | | | | | | | |
| 141 | | Analog Accel Calibration Offset Y (MSB) | | | | | | | | |
| 142 | | Analog Accel Calibration Offset Y (LSB) | | | | | | | | |
| 143 | | Analo | g Accel | Calibra | ation O | ffset Z | (MSB) | | | |
| 144 | | Analo | g Acce | l Calibr | ation C | ffset Z | (LSB) | | | |
| 145 | | Analo | g Acce | l Calibr | ation G | Gain X (| MSB) | | | |
| 146 | | Anal | og Acce | el Calib | ration (| Gain X | (LSB) | | | |
| 147 | | Analo | g Acce | l Calibr | ation G | ain Y (| MSB) | | | |
| 148 | | Anal | og Acc | el Calib | ration (| Gain Y | (LSB) | | | |
| 149 | | Analo | g Acce | l Calibr | ation 0 | Gain Z (| MSB) | | | |
| 150 | | Anal | og Acc | el Calib | ration (| Gain Z | (LSB) | | | |
| 151 | | An | alog Ad | ccel Cal | ibratio | n Align | XX | | | |
| 152 | | An | alog Ad | ccel Cal | ibratio | n Align | XY | | | |
| 153 | | An | alog A | ccel Cal | ibratio | n Align | XZ | | | |
| 154 | | Analog Accel Calibration Align YX | | | | | | | | |
| 155 | Analog Accel Calibration Align YY | | | | | | | | | |
| 156 | Analog Accel Calibration Align YZ | | | | | | | | | |
| 157 | | Analog Accel Calibration Align ZX | | | | | | | | |
| 158 | | An | alog A | ccel Cal | ibratio | n Align | ZY | | | |
| 159 | | An | alog A | ccel Cal | ibratio | n Align | ZZ | | | |



Byte # 160 - 181: Temperature (BMP180) and Pressure Calibration

| Byte | Bit 7 | Bit 6 | Bit 5 | Bit 4 | Bit 3 | Bit 2 | Bit 1 | Bit 0 | | |
|------|---------------------------------|---------------------------------|---------|----------|----------|---------|-------|-------|--|--|
| 160 | Temp & Pres Calibration AC1_MSB | | | | | | | | | |
| 161 | Temp & Pres Calibration AC1_LSB | | | | | | | | | |
| 162 | Temp & Pres Calibration AC2_MSB | | | | | | | | | |
| 163 | Temp & Pres Calibration AC2_LSB | | | | | | | | | |
| 164 | | Ten | np & Pr | es Calil | oration | AC3_N | ЛSВ | | | |
| 165 | | Ter | np & P | res Cali | bration | n AC3_I | LSB | | | |
| 166 | | Ten | np & Pr | es Calil | oration | AC4_N | ЛSВ | | | |
| 167 | | Ter | np & P | res Cali | bration | n AC4_I | LSB | | | |
| 168 | | Ten | np & Pr | es Calil | oration | AC5_N | ИSВ | | | |
| 169 | | Ter | np & P | res Cali | bration | n AC5_I | LSB | | | |
| 170 | | Ten | np & Pr | es Calil | oration | AC6_N | ИSВ | | | |
| 171 | | Temp & Pres Calibration AC6_LSB | | | | | | | | |
| 172 | | Tei | mp & P | res Cal | ibratio | n B1_N | 1SB | | | |
| 173 | | Temp & Pres Calibration B1_LSB | | | | | | | | |
| 174 | | Tei | որ & P | res Cal | ibratio | n B2_N | 1SB | | | |
| 175 | | Te | mp & F | res Ca | libratio | n B2_L | SB | | | |
| 176 | Temp & Pres Calibration MB_MSB | | | | | | | | | |
| 177 | Temp & Pres Calibration MB_LSB | | | | | | | | | |
| 178 | Temp & Pres Calibration MC_MSB | | | | | | | | | |
| 179 | | Temp & Pres Calibration MC_LSB | | | | | | | | |
| 180 | Temp & Pres Calibration MD_MSB | | | | | | | | | |
| 181 | | Ter | np & P | res Cal | ibratio | n MD_L | SB | | | |

6.3. LogAndStream source code header file (Shimmer.h)

```
/*
 * Copyright (c) 2013, Shimmer Research, Ltd.
 * All rights reserved
 *
 * Redistribution and use in source and binary forms, with or without
 * modification, are permitted provided that the following conditions are
 * met:
 *
 * Redistributions of source code must retain the above copyright
```

- * Redistributions of source code must retain the above copyright notice, this list of conditions and the following disclaimer.
- * Redistributions in binary form must reproduce the above copyright notice, this list of conditions and the following disclaimer in the documentation and/or other materials provided with the distribution.
- * Neither the name of Shimmer Research, Ltd. nor the names of its contributors may be used to endorse or promote products derived from this software without specific prior written permission.
- * You may not use or distribute this Software or any derivative works
 in any form for commercial purposes with the exception of commercial
 purposes when used in conjunction with Shimmer products purchased
 from Shimmer or their designated agent or with permission from
- * Shimmer.
- * Examples of commercial purposes would be running business
- * operations, licensing, leasing, or selling the Software, or
 - distributing the Software for use with commercial products.



```
* THIS SOFTWARE IS PROVIDED BY THE COPYRIGHT HOLDERS AND CONTRIBUTORS
 * "AS IS" AND ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING, BUT NOT
 * LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR
 * A PARTICULAR PURPOSE ARE DISCLAIMED. IN NO EVENT SHALL THE COPYRIGHT
 * OWNER OR CONTRIBUTORS BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL,
 * SPECIAL, EXEMPLARY, OR CONSEQUENTIAL DAMAGES (INCLUDING, BUT NOT
 * LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS OR SERVICES; LOSS OF USE,
 * DATA, OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER CAUSED AND ON ANY
 * THEORY OF LIABILITY, WHETHER IN CONTRACT, STRICT LIABILITY, OR TORT
 st (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF THE USE
 * OF THIS SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGE.
 */
#ifndef SHIMMER BTSD H
#define SHIMMER_BTSD_H
//these are defined in the Makefile for BtStream (TinyOS)
#define DEVICE VER
                                   //Represents SR30. 0-3 for shimmer1 to shimmer3
#define FW_IDENTIFIER
                                   //Two byte firmware identifier number: 3 for BTSD,
                                   //2 for SDLog, 1 for BTStream,
#define FW VER MAJOR
                                   //Maor version number: 0-65535
#define FW VER MINOR 0
                                   //Minor version number: 0-255
#define FW_VER_INTERNAL
                                   //internal version number: 0-255
typedef enum {
       BTSD_DOCKED = 0,
       BTSD UNDOCKED BT,
       BTSD UNDOCKED SYNC,
       BTSD_STREAMING,
       BTSD STREAMING LOGGING,
       BTSD LOGGING,
       BTSD_WAITING_USR
} BTSD STATE;
typedef uint8_t bool;
#define TRUE 1
#define FALSE 0
typedef uint8_t error_t;
#define SUCCESS 1
#define FAIL 0
// Packet Types// Packet Types
#define DATA PACKET
                                                       0x00
#define INQUIRY_COMMAND
                                                       0x01
#define INQUIRY_RESPONSE
                                                       0x02
#define GET_SAMPLING_RATE_COMMAND
                                                       0x03
#define SAMPLING_RATE_RESPONSE
                                                       0x04
                                                       0x05
#define SET_SAMPLING_RATE_COMMAND
#define TOGGLE LED COMMAND
                                                       0x06
#define START_STREAMING_COMMAND
                                                             //maintain compatibility
                                                       0x07
with Shimmer2/2r BtStream
#define SET_SENSORS_COMMAND
                                                       0x08
#define SET_LSM303DLHC_ACCEL_RANGE_COMMAND
                                                       0x09
#define LSM303DLHC_ACCEL_RANGE_RESPONSE
                                                       0x0A
#define GET_LSM303DLHC_ACCEL_RANGE_COMMAND
                                                       0x0B
#define SET_CONFIG_SETUP_BYTES_COMMAND
                                                       0x0E
#define CONFIG_SETUP_BYTES_RESPONSE
                                                       0x0F
#define GET_CONFIG_SETUP_BYTES_COMMAND
                                                       0x10
#define SET_A_ACCEL_CALIBRATION_COMMAND
                                                       0x11
#define A_ACCEL_CALIBRATION_RESPONSE
                                                       0x12
#define GET_A_ACCEL_CALIBRATION_COMMAND
                                                       0x13
```



```
#define SET_MPU9150_GYRO_CALIBRATION_COMMAND
                                                       0x14
#define MPU9150 GYRO CALIBRATION RESPONSE
                                                       0x15
#define GET_MPU9150_GYRO_CALIBRATION_COMMAND
                                                       0x16
#define SET_LSM303DLHC_MAG_CALIBRATION_COMMAND
                                                       0x17
#define LSM303DLHC_MAG_CALIBRATION_RESPONSE
                                                       0x18
#define GET_LSM303DLHC_MAG_CALIBRATION_COMMAND
                                                       0x19
#define SET_LSM303DLHC_ACCEL_CALIBRATION_COMMAND
                                                       0x1A
#define LSM303DLHC_ACCEL_CALIBRATION_RESPONSE
                                                       0x1B
#define GET_LSM303DLHC_ACCEL_CALIBRATION_COMMAND
                                                       0x1C
#define STOP_STREAMING_COMMAND
                                                       0x20//maintain compatibility
                                                        //with Shimmer2/2r BtStream
#define SET_GSR_RANGE_COMMAND
                                                       0x21
#define GSR_RANGE_RESPONSE
                                                       0x22
#define GET_GSR_RANGE_COMMAND
                                                       0x23
#define DEPRECATED_GET_DEVICE_VERSION_COMMAND
                                                       0x24
                                   //deprecated because 0x24 ('$' ASCII) as a command
                                   //is problematic if remote config is enabled in
                                   //RN42 Bluetooth module. Replaced with 0x3F command
#define DEVICE_VERSION_RESPONSE
                                                       0x25
                                   //maintain compatibility with Shimmer2/2r BtStream
#define GET ALL CALIBRATION COMMAND
                                                       0x2C
#define ALL CALIBRATION RESPONSE
                                                       0x2D
#define GET_FW_VERSION_COMMAND
                                                       0x2E
                                   //maintain compatibility with Shimmer2/2r BtStream
#define FW_VERSION_RESPONSE
                                                       0x2F
                                   //maintain compatibility with Shimmer2/2r BtStream
#define SET CHARGE STATUS LED COMMAND
                                                       0x30
#define CHARGE STATUS LED RESPONSE
                                                       0x31
#define GET_CHARGE_STATUS_LED_COMMAND
                                                       0x32
#define BUFFER_SIZE_RESPONSE
                                                       0x35
#define GET BUFFER SIZE COMMAND
                                                       0x36
#define SET LSM303DLHC MAG GAIN COMMAND
                                                       0x37
#define LSM303DLHC_MAG_GAIN_RESPONSE
                                                       0x38
#define GET_LSM303DLHC_MAG_GAIN_COMMAND
                                                       0x39
#define SET_LSM303DLHC_MAG_SAMPLING_RATE_COMMAND
                                                       0x3A
#define LSM303DLHC_MAG_SAMPLING_RATE_RESPONSE
                                                       0x3B
#define GET_LSM303DLHC_MAG_SAMPLING_RATE_COMMAND
                                                       0x3C
#define UNIQUE_SERIAL_RESPONSE
                                                       0x3D
#define GET_UNIQUE_SERIAL_COMMAND
                                                       0x3E
#define GET_DEVICE_VERSION_COMMAND
                                                       0x3F
#define SET_LSM303DLHC_ACCEL_SAMPLING_RATE_COMMAND
                                                       0x40
#define LSM303DLHC_ACCEL_SAMPLING_RATE_RESPONSE
                                                       0x41
#define GET_LSM303DLHC_ACCEL_SAMPLING_RATE_COMMAND
                                                       0x42
#define SET_LSM303DLHC_ACCEL_LPMODE_COMMAND
                                                       0x43
#define LSM303DLHC_ACCEL_LPMODE_RESPONSE
                                                       0x44
#define GET_LSM303DLHC_ACCEL_LPMODE_COMMAND
                                                       0x45
#define SET_LSM303DLHC_ACCEL_HRMODE_COMMAND
                                                       0x46
#define LSM303DLHC_ACCEL_HRMODE_RESPONSE
                                                       0x47
#define GET_LSM303DLHC_ACCEL_HRMODE_COMMAND
                                                       0x48
#define SET MPU9150 GYRO RANGE COMMAND
                                                       0x49
#define MPU9150 GYRO RANGE RESPONSE
                                                       0x4A
#define GET MPU9150 GYRO RANGE COMMAND
                                                       0x4B
#define SET_MPU9150_SAMPLING_RATE_COMMAND
                                                       0x4C
#define MPU9150_SAMPLING_RATE_RESPONSE
                                                       0x4D
#define GET_MPU9150_SAMPLING_RATE_COMMAND
                                                       0x4E
#define SET_MPU9150_ACCEL_RANGE_COMMAND
                                                       0x4F
#define MPU9150_ACCEL_RANGE_RESPONSE
                                                       0x50
#define GET_MPU9150_ACCEL_RANGE_COMMAND
                                                       0x51
#define SET_BMP180_PRES_OVERSAMPLING_RATIO_COMMAND
                                                       0x52
#define BMP180_PRES_OVERSAMPLING_RATIO_RESPONSE
                                                       0x53
#define GET_BMP180_PRES_OVERSAMPLING_RATIO_COMMAND
                                                       0x54
#define BMP180_CALIBRATION_COEFFICIENTS_RESPONSE
                                                       0x58
```



```
#define GET_BMP180_CALIBRATION_COEFFICIENTS_COMMAND
                                                        0x59
#define RESET_TO_DEFAULT_CONFIGURATION_COMMAND
                                                        0x5A
#define RESET_CALIBRATION_VALUE_COMMAND
                                                        0x5B
#define MPU9150_MAG_SENS_ADJ_VALS_RESPONSE
                                                        0x5C
#define GET_MPU9150_MAG_SENS_ADJ_VALS_COMMAND
                                                        0x5D
#define SET_INTERNAL_EXP_POWER_ENABLE_COMMAND
                                                        0x5E
#define INTERNAL_EXP_POWER_ENABLE_RESPONSE
                                                        0x5F
#define GET_INTERNAL_EXP_POWER_ENABLE_COMMAND
                                                        0x60
#define SET_EXG_REGS_COMMAND
                                                        0x61
#define EXG_REGS_RESPONSE
                                                        0x62
#define GET_EXG_REGS_COMMAND
                                                        0x63
#define SET_DAUGHTER_CARD_ID_COMMAND
                                                        0x64
#define DAUGHTER_CARD_ID_RESPONSE
                                                        0x65
#define GET_DAUGHTER_CARD_ID_COMMAND
                                                        0x66
#define SET_DAUGHTER_CARD_MEM_COMMAND
                                                        0x67
#define DAUGHTER_CARD_MEM_RESPONSE
                                                        0x68
#define GET_DAUGHTER_CARD_MEM_COMMAND
                                                        0x69
#define ACK_COMMAND_PROCESSED
                                                        0xFF
#define ROUTINE COMMUNICATION
                                                                       0xE0
                                                                              //'0'
#define START SDBT COMMAND
                                                                       0x70
#define STATUS RESPONSE
                                                                       0x71
#define GET STATUS COMMAND
                                                                       0x72
#define SET_TRIAL_CONFIG_COMMAND
                                                                       0x73
#define TRIAL_CONFIG_RESPONSE
                                                                       0x74
#define GET_TRIAL_CONFIG_COMMAND
                                                                       0x75
#define SET_CENTER_COMMAND
                                                                       0x76
#define CENTER RESPONSE
                                                                       0x77
#define GET CENTER COMMAND
                                                                       0x78
#define SET_SHIMMERNAME_COMMAND
                                                                       0x79
#define SHIMMERNAME RESPONSE
                                                                       0x7a
#define GET_SHIMMERNAME_COMMAND
                                                                       0x7b
#define SET_EXPID_COMMAND
                                                                       0x7c
#define EXPID RESPONSE
                                                                       0x7d
#define GET EXPID COMMAND
                                                                       0x7e
#define SET_MYID_COMMAND
                                                                       0x7F
#define MYID_RESPONSE
                                                                       0x80
#define GET_MYID_COMMAND
                                                                       0x81
#define SET_NSHIMMER_COMMAND
                                                                       0x82
#define NSHIMMER RESPONSE
                                                                       0x83
#define GET_NSHIMMER_COMMAND
                                                                       0x84
#define SET_CONFIGTIME_COMMAND
                                                                       0x85
#define CONFIGTIME_RESPONSE
                                                                       0x86
#define GET_CONFIGTIME_COMMAND
                                                                       0x87
#define DIR_RESPONSE
                                                                       0x88
#define GET_DIR_COMMAND
                                                                       0x89
#define INSTREAM CMD RESPONSE
                                                                       0x8A
//SENSORS0
#define SENSOR A ACCEL
                                                         0x80
#define SENSOR MPU9150 GYRO
                                                         0x40
#define SENSOR LSM303DLHC MAG
                                                         0x20
#define SENSOR EXG1 24BIT
                                                         0x10
#define SENSOR_EXG2_24BIT
                                                         0x08
#define SENSOR_GSR
                                                         0x04
#define SENSOR_EXT_A7
                                                         0x02
#define SENSOR_EXT_A6
                                                         0x01
//SENSORS1
                                                         0x80
//#define SDH_SENSOR_STRAIN
//#define SDH_SENSOR_HR
                                                         0x40
#define SENSOR_VBATT
                                                         0x20
#define SENSOR_LSM303DLHC_ACCEL
                                                         0x10
#define SENSOR_EXT_A15
                                                         0x08
```



```
#define SENSOR_INT_A1
                                                         0x04
#define SENSOR_INT_A12
                                                         0x02
#define SENSOR_INT_A13
                                                         0x01
//SENORS2
#define SENSOR_INT_A14
                                                         0x80
#define SENSOR_MPU9150_ACCEL
                                                         0x40
#define SENSOR_MPU9150_MAG
                                                         0x20
#define SENSOR_EXG1_16BIT
                                                         0x10
#define SENSOR_EXG2_16BIT
                                                         0x08
#define SENSOR_BMP180_PRESSURE
                                                         0x04
//#define SENSOR_BMP180_TEMPERATURE
                                                         0x02
#define MAX_COMMAND_ARG_SIZE
                                        //maximum number of arguments for any command
                                  131
                                        //sent
                                        //(daughter card mem write)
#define RESPONSE_PACKET_SIZE
                                  131
                                        //biggest possibly required (daughter card mem
                                        //read + 1 byte for ack)
#define MAX_NUM_CHANNELS
                                  28
                                        //3xanalogAccel + 3xdigiGyro + 3xdigiMag +
                                        //3xLSM303DLHCAccel + 3xMPU9150Accel +
                                        //3xMPU9150MAG + BMP180TEMP + BMP180PRESS +
                                        //batteryVoltage + 3xexternalADC +
                                        //4xinternalADC
#define DATA_PACKET_SIZE
                                  84
                                        //3 + (MAX_NUM_CHANNELS * 2) + 1 + 6 (+1 as
                                        //BMP180 pressure requires 3 bytes, +6 for 4
                                        //(3 byte) ExG channels plus 2 status bytes
                                        //instead of 4xinternalADC)
// Channel contents
#define X_A_ACCEL
                                           0x00
#define Y_A_ACCEL
                                           0x01
#define Z_A_ACCEL
                                           0x02
#define VBATT
                                           0x03
#define X LSM303DLHC ACCEL
                                           0x04
#define Y_LSM303DLHC_ACCEL
                                           0x05
#define Z_LSM303DLHC_ACCEL
                                           0x06
#define X_LSM303DLHC_MAG
                                           0x07
#define Y_LSM303DLHC_MAG
                                           0x08
#define Z_LSM303DLHC_MAG
                                           0x09
#define X_MPU9150_GYRO
                                           A0x0
#define Y_MPU9150_GYRO
                                           0x0B
#define Z_MPU9150_GYRO
                                           0x0C
#define EXTERNAL_ADC_7
                                           0x0D
#define EXTERNAL_ADC_
                                           0x0E
#define EXTERNAL_ADC_15
                                           0x0F
#define INTERNAL_ADC_
                                           0x10
#define INTERNAL_ADC_12
                                           0x11
#define INTERNAL_ADC_13
                                           0x12
#define INTERNAL_ADC_
                                           0x13
#define X MPU9150 ACCEL
                                           0x14
#define Y_MPU9150_ACCEL
                                           0x15
#define Z MPU9150 ACCEL
                                           0x16
#define X MPU9150 MAG
                                           0x17
#define Y_MPU9150_MAG
                                           0x18
#define Z_MPU9150_MAG
                                           0x19
#define BMP180_TEMP
                                           0x1A
#define BMP180_PRESSURE
                                           0x1B
#define GSR_RAW
                                           0x1C
#define EXG_ADS1292R_1_STATUS
                                           0x1D
#define EXG_ADS1292R_1_CH1_24BIT
                                           0x1E
#define EXG_ADS1292R_1_CH2_24BIT
                                           0x1F
#define EXG_ADS1292R_2_STATUS
                                           0x20
#define EXG_ADS1292R_2_CH1_24BIT
                                           0x21
```



```
#define EXG_ADS1292R_2_CH2_24BIT
                                           0x22
#define EXG_ADS1292R_1_CH1_16BIT
                                           0x23
#define EXG_ADS1292R_1_CH2_16BIT
                                           0x24
#define EXG_ADS1292R_2_CH1_16BIT
                                           0x25
#define EXG_ADS1292R_2_CH2_16BIT
                                           0x26
// Infomem contents
//#define NV_NUM_CONFIG_BYTES
                                          100
// Infomem contents
#define NV_NUM_SETTINGS_BYTES
                                           33
#define NV_NUM_CALIBRATION_BYTES
                                           84
#define NV_TOTAL_NUM_CONFIG_BYTES
NV_NUM_SETTINGS_BYTES+NV_NUM_CALIBRATION_BYTES
#define NV_SAMPLING_RATE
                                          0
#define NV_BUFFER_SIZE
                                          2
#define NV_SENSORS0
                                          3
#define NV_SENSORS1
                                          4
#define NV_SENSORS2
                                          5
#define NV_CONFIG_SETUP_BYTE0
                                          6
                                                 //sensors setting bytes
#define NV CONFIG SETUP BYTE1
                                          7
#define NV CONFIG SETUP BYTE2
                                          8
#define NV_CONFIG_SETUP_BYTE3
                                          9
#define NV_TRIAL_CONFIG0
                                          10
#define NV_TRIAL_CONFIG1
                                                 //debug
                                          11
#define NV_BT_INTERVAL
                                          12
#define NV EXG ADS1292R 1 CONFIG1
                                           13
                                                 // exg bytes, not implemented yet
#define NV EXG ADS1292R 1 CONFIG2
                                           14
#define NV_EXG_ADS1292R_1_LOFF
                                           15
#define NV_EXG_ADS1292R_1_CH1SET
                                           16
#define NV_EXG_ADS1292R_1_CH2SET
                                           17
#define NV_EXG_ADS1292R_1_RLD_SENS
                                           18
#define NV_EXG_ADS1292R_1_LOFF_SENS
                                           19
#define NV_EXG_ADS1292R_1_LOFF_STAT
                                           20
                                           21
#define NV_EXG_ADS1292R_1_RESP1
                                           22
#define NV_EXG_ADS1292R_1_RESP2
#define NV_EXG_ADS1292R_2_CONFIG1
                                           23
#define NV_EXG_ADS1292R_2_CONFIG2
                                           24
                                           25
#define NV_EXG_ADS1292R_2_LOFF
#define NV_EXG_ADS1292R_2_CH1SET
                                           26
#define NV_EXG_ADS1292R_2_CH2SET
                                           27
#define NV_EXG_ADS1292R_2_RLD_SENS
                                           28
#define NV_EXG_ADS1292R_2_LOFF_SENS
                                           29
#define NV_EXG_ADS1292R_2_LOFF_STAT
                                           30
#define NV_EXG_ADS1292R_2_RESP1
                                           31
#define NV_EXG_ADS1292R_2_RESP2
                                           32
#define NV_A_ACCEL_CALIBRATION
                                           33
#define NV_MPU9150_GYRO_CALIBRATION
                                           54
#define NV_LSM303DLHC_MAG_CALIBRATION
                                           75
#define NV LSM303DLHC ACCEL CALIBRATION
//Config byte masks
//Config Byte0
#define LSM303DLHC_ACCEL_SAMPLING_RATE
                                                         0xF0
#define LSM303DLHC_ACCEL_RANGE
                                                         0x0C
#define LSM303DLHC_ACCEL_LOW_POWER_MODE
                                                         0x02
#define LSM303DLHC_ACCEL_HIGH_RESOLUTION_MODE
                                                 0x01
//Config Byte1
#define MPU9150_SAMPLING_RATE
                                                         0xFF
//Config Byte2
#define LSM303DLHC_MAG_GAIN
                                                         0xE0
#define LSM303DLHC_MAG_SAMPLING_RATE
                                                         0x1C
```



```
#define MPU9150_GYRO_RANGE
                                                        0x03
//Config Byte3
#define MPU9150_ACCEL_RANGE
                                                        0xC0
#define BMP180_PRESSURE_RESOLUTION
                                                        0x30
#define GSR_RANGE
                                                        0x0E
#define EXP_POWER_ENABLE
                                                        0x01
//Unused bits 3-0
//ADC initialisation mask
#define MASK_A_ACCEL
                            0x0001
#define MASK_VBATT
                            0x0002
#define MASK_EXT_A7
                            0x0004
#define MASK_EXT_A6
                            0x0008
#define MASK_EXT_A15
                            0x0010
#define MASK_INT_A1
                            0x0020
#define MASK_INT_A12
                            0x0040
#define MASK_INT_A13
                            0x0080
#define MASK_INT_A14
                            0x0100
#define MASK_MSP_TEMP
                            0x0200
//LSM303DLHC Accel Range
//Corresponds to the FS field of the LSM303DLHC's CTRL REG4 A register
//and the AFS_SEL field of the MPU9150's ACCEL_CONFIG register
#define ACCEL_2G
                            0x00
#define ACCEL_4G
                            0x01
#define ACCEL_8G
                            0x02
#define ACCEL 16G
                            0x03
//LSM303DLHC Accel Sampling Rate
//Corresponds to the ODR field of the LSM303DLHC's CTRL REG1 A register
#define LSM303DLHC_ACCEL_POWER_DOWN
                                          0x00
#define LSM303DLHC_ACCEL_1HZ
                                          0x01
#define LSM303DLHC_ACCEL_10HZ
                                          0x02
#define LSM303DLHC_ACCEL_25HZ
                                          0x03
#define LSM303DLHC_ACCEL_50HZ
                                          0x04
#define LSM303DLHC_ACCEL_100HZ
                                          0x05
#define LSM303DLHC_ACCEL_200HZ
                                          0x06
#define LSM303DLHC_ACCEL_400HZ
                                          0x07
#define LSM303DLHC_ACCEL_1_620KHZ
                                          0x08
                                                 //1.620kHz in Low-power mode only
                                                 //1.344kHz in normal mode, 5.376kHz in
#define LSM303DLHC_ACCEL_1_344kHz
                                          0x09
                                                 //low-power mode
//LSM303DLHC Mag gain
                                          0x01 //+/-1.3 Gauss
#define LSM303DLHC_MAG_1_3G
                                          0x02 //+/-1.9 Gauss
#define LSM303DLHC_MAG_1_9G
                                          0x03 //+/-2.5 Gauss
#define LSM303DLHC_MAG_2_5G
                                          0x04 //+/-4.0 Gauss
#define LSM303DLHC_MAG_4_0G
                                          0x05 //+/-4.7 Gauss
#define LSM303DLHC_MAG_4_7G
                                          0x06 //+/-5.6 Gauss
#define LSM303DLHC_MAG_5_6G
#define LSM303DLHC_MAG_8_1G
                                          0x07 //+/-8.1 Gauss
//LSM303DLHC Mag sampling rate
#define LSM303DLHC_MAG_0_75HZ
                                          0x00 //0.75 Hz
#define LSM303DLHC_MAG_1_5HZ
                                          0x01 //1.5 Hz
#define LSM303DLHC_MAG_3HZ
                                          0x02 //3.0 Hz
#define LSM303DLHC_MAG_7_5HZ
                                          0x03 //7.5 Hz
#define LSM303DLHC_MAG_15HZ
                                          0x04 //15 Hz
#define LSM303DLHC MAG 30HZ
                                          0x05 //30 Hz
#define LSM303DLHC_MAG_75HZ
                                          0x06 //75 Hz
#define LSM303DLHC_MAG_220HZ
                                          0x07 //220 Hz
//calibration info
```



```
#define S_ACCEL
                                           0
#define S_GYRO
                                           1
                                           2
#define S_MAG
#define S_ACCEL_A
                                           3
//#define S_ECG
                                           3
                                           4
//#define S_EMG
//MPU9150 Gyro range
#define MPU9150_GYRO_250DPS
                                           0x00 //+/-250 dps
#define MPU9150_GYRO_500DPS
                                           0x01 //+/-500 dps
#define MPU9150_GYRO_1000DPS
                                           0x02 //+/-1000 dps
#define MPU9150_GYRO_2000DPS
                                           0x03 //+/-2000 dps
//#digital accel_range
#define RANGE_2G
                                           0
#define RANGE_4G
                                           1
#define RANGE_8G
                                           2
#define RANGE_16G
                                           3
//#mag_gain
#define LSM303 MAG 13GA
                                           1
#define LSM303 MAG 19GA
                                           2
#define LSM303_MAG_25GA
                                           3
#define LSM303_MAG_40GA
                                           4
#define LSM303_MAG_47GA
                                           5
#define LSM303_MAG_56GA
                                           6
#define LSM303 MAG 81GA
                                           7
//SD Log file header format
#define SDHEAD_LEN
                                           256// 0-255
#define SDH_SAMPLE_RATE_0
                                                  0
#define SDH_SAMPLE_RATE_1
                                                  1
                                                  3
#define SDH_SENSORS0
                                                  4
#define SDH_SENSORS1
                                                  5
#define SDH_SENSORS2
#define SDH_CONFIG_SETUP_BYTE0
                                                  8 //sensors setting bytes
                                                  9
#define SDH_CONFIG_SETUP_BYTE1
#define SDH_CONFIG_SETUP_BYTE2
                                                  10
#define SDH_CONFIG_SETUP_BYTE3
                                                  11
#define SDH_TRIAL_CONFIG0
                                                  16
#define SDH_TRIAL_CONFIG1
                                                  17
#define SDH_BROADCAST_INTERVAL
                                                  18
                                                  22
//#define SDH_ACCEL_RANGE
                                                         //digital accel range
                                                  23
//#define SDH_GSR_RANGE
//#define SDH_MAG_RANGE
                                                  24
//#define SDH_MAG_RATE
                                                  25
//#define SDH_ACCEL_RATE
                                                  26
                                                         digital accel rate
//#define SDH_GYRO_RATE
                                                  27
//#define SDH_GYRO_RANGE
                                                  28
//#define SDH PRESSURE PREC
                                                  29
#define SDH SHIMMERVERSION BYTE 0
                                                  30
#define SDH_SHIMMERVERSION_BYTE_1
                                                  31
#define SDH_MYTRIAL_ID
                                                  32
#define SDH_NSHIMMER
                                                  33
#define SDH_FW_VERSION_TYPE_0
                                                  34
#define SDH_FW_VERSION_TYPE_1
                                                  35
#define SDH_FW_VERSION_MAJOR_0
                                                  36
#define SDH_FW_VERSION_MAJOR_1
                                                  37
#define SDH_FW_VERSION_MINOR
                                                  38
#define SDH_FW_VERSION_INTERNAL
                                                  39
#define SDH_CONFIG_TIME_0
                                                  52
```



```
#define SDH_CONFIG_TIME_1
                                                  53
#define SDH_CONFIG_TIME_
                                                  54
                                                  55
#define SDH_CONFIG_TIME_3
#define SDH_EXG_ADS1292R_1_CONFIG1
                                                  56
#define SDH_EXG_ADS1292R_1_CONFIG2
                                                  57
#define SDH_EXG_ADS1292R_1_LOFF
                                                  58
#define SDH_EXG_ADS1292R_1_CH1SET
                                                  59
#define SDH_EXG_ADS1292R_1_CH2SET
                                                  60
#define SDH_EXG_ADS1292R_1_RLD_SENS
                                                  61
#define SDH_EXG_ADS1292R_1_LOFF_SENS
                                                  62
#define SDH_EXG_ADS1292R_1_LOFF_STAT
                                                  63
#define SDH_EXG_ADS1292R_1_RESP1
                                                  64
#define SDH_EXG_ADS1292R_1_RESP2
                                                  65
#define SDH_EXG_ADS1292R_2_CONFIG1
                                                  66
#define SDH_EXG_ADS1292R_2_CONFIG2
                                                  67
#define SDH_EXG_ADS1292R_2_LOFF
                                                  68
#define SDH_EXG_ADS1292R_2_CH1SET
                                                  69
#define SDH_EXG_ADS1292R_2_CH2SET
                                                  70
#define SDH_EXG_ADS1292R_2_RLD_SENS
                                                  71
#define SDH_EXG_ADS1292R_2_LOFF_SENS
                                                  72
                                                  73
#define SDH_EXG_ADS1292R_2_LOFF_STAT
                                                  74
#define SDH EXG ADS1292R 2 RESP1
#define SDH_EXG_ADS1292R_2_RESP2
                                                  75
#define SDH_LSM303DLHC_ACCEL_CALIBRATION
                                                  76
#define SDH_MPU9150_GYRO_CALIBRATION
                                                  97
#define SDH_LSM303DLHC_MAG_CALIBRATION
                                                  118
#define SDH A ACCEL CALIBRATION
                                                  139
#define SDH TEMP PRES CALIBRATION
                                                  160
#define SDH_MY_LOCALTIME
                                                  252
                                                         //252-255
//SENSORS0
#define
              SDH SENSOR A ACCEL
                                                  0x80
#define SDH_SENSOR_MPU9150_GYRO
                                                  0x40
#define SDH_SENSOR_LSM303DLHC_MAG
                                                  0x20
#define SDH_SENSOR_EXG1_24BIT
                                                  0x10
#define SDH_SENSOR_EXG2_24BIT
                                                  0x08
#define SDH_SENSOR_GSR
                                                  0x04
#define SDH_SENSOR_EXTCH7
                                                  0x02
#define SDH_SENSOR_EXTCH6
                                                  0x01
//SENSORS1
//#define SDH_SENSOR_STRAIN
                                                  0x80
//#define SDH_SENSOR_HR
                                                  0x40
#define SDH_SENSOR_VBATT
                                                  0x20
#define SDH_SENSOR_LSM303DLHC_ACCEL
                                                  0x10
#define SDH_SENSOR_EXTCH15
                                                  0x08
#define SDH_SENSOR_INTCH1
                                                  0x04
#define SDH_SENSOR_INTCH12
                                                  0x02
#define SDH_SENSOR_INTCH13
                                                  0x01
//SENSORS2
#define SDH SENSOR INTCH14
                                                  0x80
#define SDH SENSOR MPU9150 ACCEL
                                                  0x40
#define SDH SENSOR MPU9150 MAG
                                                  0x20
#define SDH_SENSOR_EXG1_16BIT
                                                  0x10
#define SDH_SENSOR_EXG2_16BIT
                                                  0x08
#define SDH_SENSOR_BMP180_PRES
                                                  0x04
//#define SDH_SENSOR_BMP180_TEMP
                                                  0x02
//SENSORS3
#define SDH_SENSOR_MSP430_TEMP
                                                  0x01
#define SDH_SENSOR_TCXO
                                                  0x80
//SDH_TRIAL_CONFIG0
#define SDH_IAMMASTER
                                                  0x02
```



```
#define SDH_TIME_SYNC
                                                 0x04
#define SDH_TIME_STAMP
                                                 0x08// not used now, reserved as 1
#define SDH_GYRO_BUTTON_ENABLE
#define SDH_USER_BUTTON_ENABLE
                                                 0x20
#define SDH_SET_PMUX
                                                 0x40// not used now, reserved as 0
#define SDH_SET_5V_REG
                                                 0x80// not used now
//SDH_TRIAL_CONFIG1
#define SDH_SINGLETOUCH
                                                 0x80
#define SDH_ACCEL_LPM
                                                 0x40//config has this bit
#define SDH_ACCEL_HRM
                                                 0x20//config has this bit
#define SDH_TCXO
#define SDH_EXP_POWER
                                                 0x08//config has this bit
#define SDH_MONITOR
                                                 0x04
//choice of clock
#define TCXO CLOCK
                           255765.625
#define MSP430_CLOCK 32768
// BT routine communication
// all node time must *2 in use
// all center time must *4 in use
#define RC AHD
#define RC_WINDOW_N
                           13
#define RC_WINDOW_C
                           27
#define RC_INT_N
                          27
#define RC_INT_C
                           54
                                         //240
#define RC CLK N
                           16384
                                         //16384=2hz;//32768=1hz;8192=4hz
#define RC CLK C
                           8192
                                         //16384=2hz;//32768=1hz;8192=4hz
#define RC_FACTOR_N 32768/RC_CLK_N
                                         //16384=2hz;//32768=1hz;8192=4hz
#define RC_FACTOR_C 32768/RC_CLK_C
                                         //16384=2hz;//32768=1hz;8192=4hz
//routine communication response text - ack:flag:time4:time3:time2:time1
#define RCT SIZE
                           6
                           0
#define RCT ACK
                           1
#define RCT_FLG
                           2
#define RCT_TIME
#define RCT_TIME_SIZE
// sd card write buffer size
                                  //4095 255
#define SDBUFF_SIZE_MAX
                           4096
#define SDBUFF_SIZE
                           512
                                   //4095 512
// BATTERY
                           0
#define BATT HIGH
#define BATT_MID
                           1
#define BATT_LOW
                           19660800
#define BATT_ITNERVAL
#define MAX_CHARS
//BMP Pressure oversampling ratio
#define BMP180 OSS 1
                                   0x00
#define BMP180_OSS_2
                                   0x01
#define BMP180_OSS_4
                                   0x02
#define BMP180_OSS_8
                                   0x03
//BtStream specific extension to range values : should SDLog keep it?
#define GSR_AUTORANGE
                                   0x04
#endif
```



6.4. Troubleshoot

Green and Blue LEDs flash when I undock or reboot my Shimmer.

The green and blue LEDs flashing on the *Shimmer3* indicate that the firmware either cannot create the required directories or cannot write to the required files on the SD card. Try the following steps to rectify the problem:

- 1. Check that the SD card is correctly inserted.
- 2. Ensure that the SD card has a capacity of 2 GB or less and is compatible with *Shimmer3* (refer to the *Shimmer User Manual*, available for download from www.shimmersensing.com).
- 3. Ensure that the SD card memory is not full.
- 4. Ensure that the *experimentid* and *shimmername* parameters, specified in the *sdlog.cfg* file contain only alphanumeric characters (a,..., z, A,..., Z, 0,..., 9), dash ('-') and underscore ('_').

When I undock or power on my Shimmer, the green LED continuously flashes at a rate of 5Hz

Try power-cycling the *Shimmer3*. The problem may be due to an error in bluetooth initialisation.

If the problem persists after power-cycling, try changing the SD card for a newer one - if the SD card is corrupt, the firmware may fail to correctly read the configuration file.

The data file is empty after logging

Ensure that you log for a minimum of one minute in order for data to be written to the SD card.

The configuration does not match the parameters in the sdlog.cfg file.

The *sdlog.cfg* file is read once each time the *Shimmer3* is rebooted or undocked so the configuration will always match the most recent configuration file at the time of logging. If the parameters in the configuration header of your data files do not match those in the *sdlog.cfg* file, it is likely that you have changed the *sdlog.cfg* file contents since logging the data in question.

The calibration parameters in the configuration header do not match the Calibration/calibParams.ini file.

Calibration parameters will only be loaded for sensors that are enabled; calibration parameters for disabled sensors will all have zero value.

If a sensor is enabled and the calibration parameters do not match the calibration file, try the following steps to rectify the problem:



- 1. Ensure that the calibration file is stored in the correct file location from the SD card top level directory, according to the following (case-sensitive) options:
 - /Calibration/calibParams.ini
 - /calibration/calibParams.ini

No other file path will be recognised by the firmware.

2. Ensure that you have implemented the correct byte-order and endianness when you read the calibration parameters from the configuration header, according to the information in Section 5.2.

Shimmer

The Realtime Building Clonshaugh Technology Park

T: +353 (1)848 6 112

E: info@ShimmerSensing.com

Shimmer North America

USA **T:** +1 857 362 X254

E: info@ShimmerSensing.com

www.ShimmerSensing.com

/ShimmerResearch

in /company/Shimmer

→ ShimmerResearch

→ Shimme